


DERMABOND®
TOPICAL SKIN ADHESIVE

STRONG CLOSURE. STRONG PROTECTION.

DERMABOND® ProPen Topical Skin Adhesive 2 - Octyl Cyanoacrylate

DERMABOND® Adhesive provides microbial barrier protection as well as final skin closure in a variety of surgical procedures. You can use it with confidence on small incisions as well as on longer incisions that are easily approximated.



Protects the Patient and Enhances Patient Comfort

- It provides an effective seal that sloughs off naturally over the course of 5 to 10 days.
- Additional dressings can be used, but are not required¹.
- Patients can shower immediately¹.

Proven Protection

Studies have been conducted to evaluate the effectiveness of DERMABOND® Topical Skin Adhesive in reducing the incidence of Surgical Site Infection (SSI).

- In an in vitro model¹, data indicated that DERMABOND® Topical Skin Adhesive is an effective barrier to microbial penetration, with 95% confidence of 99% efficacy for 72 hours⁴.
- In a clinical trial² involving plastic surgery operations, wounds protected with DERMABOND® Topical Skin Adhesive showed statistically significantly fewer incidences of SSI³.

1. DERMABOND® Adhesive package insert. Somerville, NJ. ETHICON, INC., 2003.

2. Quinn J, Sutcliffe T, et al. A randomised trial comparing octylcyanoacrylate tissue adhesive and sutures in the management of lacerations. JAMA 1997; 277: 1527-1530.

3. Blondeel PNV, et al. Closure of long incisions with a new formulation of octylcyanoacrylate tissue adhesive vs commercially available methods. AJS 2004, 188: 307-313.

4. Bhende S, Rothenburger S, Spangler DJ, Dito M. In vitro assessment of microbial barrier properties of DERMABOND® Topical Skin Adhesive. Surg Infect. 2002;3(3):251-257.

5. Silvestri A, et al. Octyl-2-Cyanoacrylate Adhesive for Skin Closure and Prevention of Infection in Plastic Surgery. Aesthetic Plastic Surgery. 2006;30:659-699.



STRONG CLOSURE. STRONG PROTECTION.

The Range at a Glance

DERMABOND® Adhesive Mini	Small vial for use in the emergency room
DERMABOND® Adhesive Regular	Regular vial for use in emergency room and operating room
DERMABOND® Adhesive ProPen	Innovative presentation facilitates ease of use and improved control over adhesive placement
DERMABOND® Adhesive ProPen XL	Delivers twice the the volume of adhesive, ideal for longer incisions



Product Name	Size	Packaging	Product Code
DERMABOND® Mini Topical Skin Adhesive 2-Octyl Cyanoacrylate	0.38 mL	12 Ampoules Per Box	AHVM12
DERMABOND® Topical Skin Adhesive 2-Octyl Cyanoacrylate	0.5 mL	6 Ampoules Per Box	AHV6
DERMABOND® Topical Skin Adhesive 2-Octyl Cyanoacrylate	0.5 mL	12 Ampoules Per Box	AHV12
DERMABOND® ProPen Topical Skin Adhesive 2-Octyl Cyanoacrylate	0.5 mL	6 Ampoules Per Box	APP6
DERMABOND® ProPen XL Topical Skin Adhesive 2-Octyl Cyanoacrylate	0.75 mL	6 Ampoules Per Box	APPXL6