


# Product Data Sheet

Disposable Goggle Foam 1105

REF 23-254

Version: 29.04.2013, rev.1

Tryckt: 10.6.2013

REF	Storlek	
23-254		
CE0434		
Producerad i enlighet med legala föreskrifter för medicintekniska produkter		
   		
<b>Material:</b>		
<b>REACH föreskrifter:</b>		
<b>Egenskaper:</b>		<p>Product name: Disposable goggle foam 1105 Part number 8011527 Indications for use Disposable goggle foam for VNG systems VO425, VN415 and VF405 including b+m variants Application Use only for VNG goggles Dimension 207x131x12 mm Weight 5 g Box Set of 24 pcs Materials:Material Polyurethane foam on polyether basis. (Velour on side for fixation) Colour: Anthracite Latex free YES PVC free YES Certificate Öko-Tex Standard 100 class 1 Environment Temperature limitations 5-35°C (Operation) Store cool and dry Keep dry, Keep away from sunlight Non-sterile Recommended directions 1. Remove old foam 2. Ensure Velcro on goggle is intact 3. Carefully place new foam with the dark velour side toward the Velcro. Cautions and precautions Single patient use only US: RX only Use only with VNG systems VO425, VN415 and VF405</p>
<b>Steril produkt:</b>		Nej
<b>Hållbarhetstid:</b>		
<b>MBA/st:</b>		1
<b>Avdelningsförpackning/st:</b>		24
<b>Transportlåda/st:</b>		24