

VICRYL*

(Polyglactin 910) Braided and Monofilament

SUTURE

| | | |
|---------------------|-------------------|---|
| Description: | Braided | For easy handling and secure knot tying. |
| | Coated | For smooth passage through tissue and easy knot tie down. |
| | Synthetic | For minimal tissue reaction. |
| | Absorbable | By hydrolysis. Essentially complete in 56-70 days. |

| | | |
|--------------------------|--------------------------|--|
| Tensile strength: | Post implantation | Approximate % original strength remaining |
| | 14 days | 75% |
| | 21 days | 50% (6-0 and larger) 40% (7-0 and smaller) |
| | 28 days | 25% (6-0 and larger) |

Color: Violet or undyed.

Range: 10-0 to 2 (USP). Supplied as needed sutures and ligatures. Sizes 10-0 and 9-0 are monofilament in structure.












Indications: VICRYL* sutures are intended for use in soft tissue approximation and/or ligation, including use in ophthalmic surgery, peripheral nerve adaptation and microsurgery for vessels less than 2 mm in diameter.

Typical areas of use include:

- Fascia closure - due to the unique 28-day claim
- Subcutaneous fat
- Joint capsule
- Uterus

Full details in the Instruction for Use included in every package.








| 1/2 Circle Taper Point Needles | | | | | | | |
|--|-------|-------|-------------|-------------------|----------------|---------|--|
| Needle | Gauge | Color | Length (cm) | Other information | Number per box | Code | |
|  TF-6 6,5 mm | 8-0 | ● | 30 | | 12 | V1027G | |
| | 7-0 | ● | 30 | | 36 | V1028H | |
|  TF-1 10 mm | 5-0 | ● | 70 | | 36 | V133H | |
|  RB-2 13 mm | 6-0 | ● | 70 | | 36 | PZ2571H | |
|  TF 13 mm | 6-0 | ● | 70 | | 36 | V926H | |
| | 5-0 | ● | 70 | | 36 | V925H | |
| | 4-0 | ● | 70 | | 36 | V924H | |
|  RB-1 17 mm | 6-0 | ● | 70 | | 36 | V302H | |
| | 5-0 | ● | 70 | | 36 | V303H | |
| | 5-0 | ○ | 70 | | 36 | V213H | |
| | 4-0 | ● | 70 | | 36 | V304H | |
| | 4-0 | ○ | 70 | | 36 | V214H | |
| | 4-0 | ● | 8x45 | CONTROL RELEASE* | 24 | V734E | |
| | 4-0 | ● | 4x45 | CONTROL RELEASE* | 12 | V735G | |
| | 4-0 | ● | 8x70 | CONTROL RELEASE* | 24 | V799E | |
| | 3-0 | ● | 70 | | 36 | V305H | |
| | 3-0 | ○ | 70 | | 36 | V215H | |
| | 3-0 | ● | 8x45 | CONTROL RELEASE* | 12 | J713D | |
| 2-0 | ● | 70 | | 36 | V306H | | |
|  RB-1 RB-1 17 mm | 3-0 | ● | 70 | | 36 | V667H | |
|  JRB-1  VISI-BLACK 17 mm | 5-0 | ● | 70 | VISI-BLACK* | 36 | V3030H | |
| | 4-0 | ● | 70 | VISI-BLACK* | 36 | V3040H | |
|  SH-1 22 mm | 4-0 | ● | 70 | | 36 | V310H | |
| | 4-0 | ○ | 70 | | 36 | V426H | |
| | 3-0 | ● | 70 | | 36 | V311H | |
| | 3-0 | ○ | 70 | | 36 | V219H | |
| | 3-0 | ● | 4x45 | CONTROL RELEASE* | 36 | V782H | |
| | 2-0 | ● | 70 | | 36 | V312H | |
| | 2-0 | ○ | 70 | | 36 | V220H | |
| | 2-0 | ● | 8x45 | CONTROL RELEASE* | 24 | V792E | |
|  JB-1  VISI-BLACK 22 mm | 4-0 | ● | 70 | VISI-BLACK* | 36 | V3100H | |
| | 3-0 | ● | 70 | VISI-BLACK* | 36 | V3110H | |
| | 2-0 | ● | 70 | VISI-BLACK* | 36 | V3120H | |





VICRYL*
(Polyglactin 910) Braided and Monofilament


|  1/2 Circle Taper Point Needles | | | | | | |
|--|-------|-------|-------------|-------------------|----------------|---------|
| Needle | Gauge | Color | Length (cm) | Other information | Number per box | Code |
|  SH 26 mm | 4-0 | ● | 70 | | 36 | V315H |
| | 4-0 | ● | 8x70 | CONTROL RELEASE* | 12 | V783G |
| | 3-0 | ● | 70 | | 36 | V316H |
| | 3-0 | ○ | 70 | | 36 | V416H |
| | 3-0 | ● | 90 | | 36 | V242H |
| | 3-0 | ● | 4x45 | CONTROL RELEASE* | 36 | V780H |
| | 3-0 | ● | 8x45 | CONTROL RELEASE* | 12 | J774D |
| | 3-0 | ● | 8x70 | CONTROL RELEASE* | 12 | J784D |
| | 2-0 | ● | 70 | | 36 | V317H |
| | 2-0 | ○ | 70 | | 36 | V224H |
| | 2-0 | ● | 90 | | 36 | V243H |
| | 2-0 | ● | 8x45 | CONTROL RELEASE* | 12 | J775D |
| | 0 | ● | 70 | | 36 | V318H |
| | 0 | ○ | 70 | | 36 | V225H |
|  JB VISI-BLACK 26 mm | 4-0 | ● | 70 | VISI-BLACK* | 36 | V3150H |
| | 3-0 | ● | 70 | VISI-BLACK* | 36 | V3160H |
| | 3-0 | ○ | 70 | VISI-BLACK* | 36 | V4416H |
| | 2-0 | ● | 70 | VISI-BLACK* | 36 | V3170H |
|  MH-1 31 mm | 3-0 | ● | 70 | | 36 | V245H |
| | 2-0 | ● | 70 | | 36 | V246H |
| | 0 | ● | 70 | | 36 | V247H |
|  MH 36 mm | 4-0 | ● | 90 | | 36 | PG2592H |
| | 3-0 | ● | 70 | | 36 | V322H |
| | 3-0 | ● | 90 | | 36 | V100H |
| | 2-0 | ● | 70 | | 36 | V323H |
| | 2-0 | ○ | 70 | | 36 | V226H |
| | 2-0 | ● | 90 | | 36 | V101H |
| | 0 | ● | 70 | | 36 | V324H |
| | 0 | ● | 90 | | 36 | V102H |
| | 1 | ● | 70 | | 36 | V325H |
| | 1 | ○ | 70 | | 36 | V228H |
|  LH 40 mm | 2-0 | ● | 70 | | 36 | PH2574H |





|  1/2 Circle Taper Point Needles (Heavier Wire Diameter) | | | | | | |
|--|-------|-------|-------------|-------------------|----------------|-----------|
| Needle | Gauge | Color | Length (cm) | Other information | Number per box | Code |
|  UCL 13 mm | 3-0 | ● | 70 | | 36 | V113H |
| | 0 | ● | 70 | | 36 | V115H |
|  UCLX 16 mm | 3-0 | ● | 70 | | 36 | V13H |
| | 2-0 | ● | 70 | | 36 | V14H |
| | 0 | ● | 70 | | 36 | V15H |
| | 1 | ● | 70 | | 36 | V16H |
|  CT-3 22 mm | 3-0 | ● | 70 | | 36 | V327H |
| | 2-0 | ● | 70 | | 36 | V328H |
| | 0 | ● | 70 | | 36 | V329H |
|  CT-2 26 mm | 3-0 | ● | 70 | | 36 | V332H |
| | 2-0 | ● | 70 | | 36 | V333H |
| | 2-0 | ○ | 70 | | 36 | V289H |
| | 2-0 | ● | 4x45 | CONTROL RELEASE* | 24 | V7842E |
| | 0 | ● | 70 | | 36 | V334H |
| | 0 | ○ | 70 | | 36 | V270H |
| | 0 | ● | 8x45 | CONTROL RELEASE* | 24 | V744E |
| | 0 | ● | 4x70 | CONTROL RELEASE* | 24 | 4PK2575YE |
| | 1 | ● | 70 | | 36 | V335H |
| | 2 | ● | 70 | | 36 | V336H |
|  CT-1 36 mm | 3-0 | ● | 70 | | 36 | V338H |
| | 3-0 | ○ | 90 | | 36 | V944H |
| | 2-0 | ● | 70 | | 36 | V339H |
| | 2-0 | ● | 90 | | 36 | V345H |
| | 2-0 | ○ | 90 | | 36 | V945H |
| | 2-0 | ○ | 8x70 | CONTROL RELEASE* | 12 | J742D |
| | 0 | ● | 70 | | 36 | V340H |
| | 0 | ● | 90 | | 36 | V346H |
| | 0 | ○ | 90 | | 36 | V946H |
| | 0 | ● | 8x45 | CONTROL RELEASE* | 12 | J740D |
| | 1 | ● | 70 | | 36 | V341H |
| | 1 | ● | 90 | | 36 | V347H |
| | 1 | ● | 4x70 | CONTROL RELEASE* | 12 | V736G |
| | 2 | ● | 90 | | 36 | V348H |
| 2 | ○ | 90 | | 36 | V948H | |
|  MO-4 36 mm | 0 | ● | 70 | | 36 | V436H |
| | 1 | ● | 90 | | 36 | V255H |

VICRYL*
(Polyglactin 910) Braided and Monofilament



|  1/2 Circle Taper Point Needles (Heavier Wire Diameter) | | | | | | |
|--|--|-------|-------------|-------------------|----------------|-----------|
| Needle | Gauge | Color | Length (cm) | Other information | Number per box | Code |
|  CT 40 mm | 3-0 | ● | 70 | | 36 | V350H |
| | 2-0 | ● | 70 | | 36 | V351H |
| | 2-0 | ● | 8x45 | CONTROL RELEASE* | 12 | 8PM2544YG |
| | 0 | ● | 70 | | 36 | V352H |
| | 0 | ● | 90 | | 36 | V358H |
| | 0 | ● | 150 | LOOP | 24 | V1407E |
| | 1 | ● | 70 | | 36 | V353H |
| | 1 | ● | 90 | | 36 | V359H |
| | 1 | ● | 8x45 | CONTROL RELEASE* | 12 | J753D |
| | 2 | ● | 90 | | 36 | V360H |
| | 2 | ○ | 90 | | 36 | V960H |
| |  CTX 48 mm | 3-0 | ● | 70 | | 36 |
| 2-0 | | ● | 70 | | 36 | V363H |
| 0 | | ● | 70 | | 36 | V364H |
| 0 | | ● | 90 | | 36 | V370H |
| 0 | | ○ | 90 | | 36 | V170H |
| 0 | | ● | 8x45 | CONTROL RELEASE* | 12 | J764D |
| 1 | | ● | 70 | | 36 | V365H |
| 1 | | ● | 90 | | 36 | V371H |
| 1 | | ○ | 90 | | 36 | V366H |
| 1 | | ● | 8x45 | CONTROL RELEASE* | 12 | J765D |
| 2 | | ● | 90 | | 36 | V372H |
| 2 | | ○ | 90 | | 36 | V367H |
|  TP-1 65 mm | | 2 | ● | 4x70 | | 6 |
|  XLH 70 mm | 2-0 | ● | 70 | | 12 | V581G |
| | 0 | ● | 70 | | 12 | V582G |
| | 1 | ● | 70 | | 12 | V583G |


| 1/2 Circle ETHIGUARD® Blunt Point Needles | | | | | | |
|---|-------|-------|-------------|-------------------|----------------|-------|
| Needle | Gauge | Color | Length (cm) | Other information | Number per box | Code |
|  MHB-1 31 mm | 3-0 | ● | 75 | | 12 | W9983 |
| | 2-0 | ● | 75 | | 12 | W9984 |
|  CTB-1 36 mm | 2-0 | ● | 75 | | 12 | W9985 |
| | 2-0 | ● | 75 | Heavy Needle | 12 | W9993 |
| | 0 | ● | 75 | | 12 | W9986 |
| | 0 | ● | 90 | Heavy Needle | 12 | W9994 |
| | 1 | ● | 90 | Heavy Needle | 12 | W9995 |
|  CTB 40 mm | 0 | ● | 75 | | 12 | W9988 |
| | 0 | ● | 90 | Heavy Needle | 12 | W9996 |
| | 1 | ● | 90 | Heavy Needle | 12 | W9997 |
|  CTXB 48 mm | 0 | ● | 90 | | 12 | W9990 |
| | 1 | ● | 90 | | 12 | W9991 |


| 3/8 Circle Blunt Point Needles | | | | | | |
|--|-------|-------|-------------|-------------------|----------------|-------|
| Needle | Gauge | Color | Length (cm) | Other information | Number per box | Code |
|  BP-5 65 mm | 1 | ● | 100 | | 12 | W9391 |

| 5/8 Circle Taper Point Needles | | | | | | |
|--|-------|-------|-------------|-------------------|----------------|----------|
| Needle | Gauge | Color | Length (cm) | Other information | Number per box | Code |
|  UR-6 26 mm | 3-0 | ● | 70 | | 36 | V601H |
| | 2-0 | ● | 70 | | 36 | V602H |
| | 0 | ● | 70 | | 36 | V603H |
| | 1 | ● | 70 | | 36 | VC2575H |
|  UR-6 26 mm | 2-0 | ● | 70 | | 36 | VCC2573H |
|  UR-5 36 mm | 2-0 | ● | 70 | | 36 | V375H |
| | 0 | ● | 70 | | 36 | V376H |
| | 1 | ● | 70 | | 36 | V377H |
|  UR-4 40 mm | 2-0 | ● | 70 | | 36 | V380H |
| | 0 | ● | 70 | | 36 | V381H |

| ▼ 1/2 Circle TAPERCUT* Needles | | | | | | | |
|---|-------|-------|-------------|-------------------|----------------|------|----------|
| Needle | Gauge | Color | Length (cm) | Other information | Number per box | Code | |
|  V-37 | 40 mm | 2-0 | ● | 90 | | 36 | VS4020H |
| | | 0 | ● | 90 | | 36 | VS4000H |
|  V-37 | 40 mm | 1 | ● | 70 | | 36 | ODD2574H |
|  V-40 | 48 mm | 2-0 | ● | 90 | | 36 | VS5020H |
| | | 0 | ● | 70 | | 36 | VS5000H |
| | | 1 | ● | 150 | LOOP | 24 | V829E |
|  BPT-1 | 48 mm | 6 | ● | 75 | | 12 | V140G |
| | | 6 | ● | 150 | | 12 | OH2528G |



| ▼ 3/8 Circle TAPERCUT* Needles | | | | | | | |
|---|-------|-------|-------------|-------------------|----------------|------|-------|
| Needle | Gauge | Color | Length (cm) | Other information | Number per box | Code | |
|  V-18 | 13 mm | 6-0 | ● | 70 | | 36 | V991H |
| | | 5-0 | ● | 70 | | 36 | V990H |
|  V-4 | 17 mm | 4-0 | ○ | 70 | | 12 | V240G |
| | | 4-0 | ○ | 70 | | 36 | V240H |




| ▼ 5/8 Circle TAPERCUT* Needles | | | | | | | |
|---|-------|-------|-------------|-------------------|----------------|------|-------|
| Needle | Gauge | Color | Length (cm) | Other information | Number per box | Code | |
|  UW-10 | 10 mm | 6-0 | ● | 30 | | 12 | W9172 |




| ▼ Compound Curve TAPERCUT* Needles | | | | | | | |
|--|-------|-------|-------------|-------------------|----------------|------|--------|
| Needle | Gauge | Color | Length (cm) | Other information | Number per box | Code | |
|  DA VISI-BLACK | 25 mm | 4-0 | ● | 45 | VISI-BLACK* | 12 | 18531G |

VICRYL*

(Polyglactin 910) Braided and Monofilament

| ● Fish Hook Needles | | | | | | |
|---|-------|-------|-------------|-------------------|----------------|---------|
| Needle | Gauge | Color | Length (cm) | Other information | Number per box | Code |
|  PB-30 30 mm | 2-0 | ● | 70 | Taper point | 36 | DF2570H |
| | 1 | ● | 70 | Taper point | 36 | DF2573H |
|  V-38 31 mm | 1 | ● | 75 | TAPERCUT* | 12 | W9335 |









| ● Straight Taper Point Needles | | | | | | |
|--|-------|-------|-------------|-------------------|----------------|---------|
| Needle | Gauge | Color | Length (cm) | Other information | Number per box | Code |
|  ST-4 19 mm | 4-0 | ● | 70 | | 36 | BA2573H |
| | 3-0 | ● | 70 | | 36 | BA2574H |
|  ST-4 19 mm | 4-0 | ● | 35 | | 36 | V620H |
|  ST 48 mm | 3-0 | ● | 70 | | 36 | BC2573H |



| ▼ 1/2 Circle TAPERCUT* Needles | | | | | | |
|--|-------|-------|-------------|-------------------|----------------|-----------|
| Needle | Gauge | Color | Length (cm) | Other information | Number per box | Code |
|  V-5 17 mm | 4-0 | ● | 70 | | 12 | V994G |
| | 4-0 | ● | 70 | | 36 | V994H |
| | 4-0 | ○ | 70 | | 36 | V230H |
| | 3-0 | ● | 70 | | 36 | V276H |
|  V-7 26 mm | 3-0 | ● | 70 | | 12 | V998G |
| | 3-0 | ● | 70 | | 36 | V998H |
| | 2-0 | ● | 70 | | 36 | V999H |
| | 2-0 | ● | 6x45 | CONTROL RELEASE* | 24 | 6OB2545YE |
| | 0 | ● | 70 | | 36 | V279H |
|  V-34 36 mm | 3-0 | ● | 90 | | 36 | V516H |
| | 2-0 | ● | 90 | | 36 | V517H |
| | 0 | ● | 90 | | 36 | V518H |
| | 1 | ● | 90 | | 36 | V519H |
| | 1 | ○ | 90 | | 36 | V919H |

VICRYL*
(Polyglactin 910) Braided and Monofilament

| ▼ 1/2 Circle Reverse Cutting Needles (Heavier Wire Diameter) | | | | | | | |
|---|-------|-------|-------------|-------------------|-----------------|------|---------|
| Needle | Gauge | Color | Length (cm) | Other information | Number per box | Code | |
|  X-1 | 22 mm | 3-0 | ● | 70 | | 36 | V460H |
| | | 3-0 | ○ | 70 | | 36 | ME2673H |
| | | 2-0 | ● | 70 | | 36 | V461H |
| | | 2-0 | ○ | 75 | | 36 | JV459 |
|  OS-4 | 22 mm | 0 | ○ | 70 | HEAVY NEEDLE | 36 | V694H |
| | | 2 | ● | 75 | HEAVY NEEDLE | 12 | W9287 |
|  CP-2 | 26 mm | 3-0 | ● | 70 | | 36 | MF2572H |
| | | 2-0 | ○ | 70 | | 36 | V869H |
| | | 0 | ● | 70 | | 36 | MF2574H |
| | | 0 | ● | 4x45 | CONTROL RELEASE | 24 | V870E |
| | | 1 | ○ | 70 | | 36 | MF2675H |
|  CP-1 | 36 mm | 3-0 | ● | 70 | | 36 | V465H |
| | | 2-0 | ● | 70 | | 36 | V466H |
| | | 2-0 | ● | 90 | | 36 | V472H |
| | | 2-0 | ○ | 90 | | 36 | JV1472 |
| | | 0 | ● | 70 | | 36 | V467H |
| | | 0 | ● | 90 | | 36 | V473H |
| | | 0 | ● | 4x45 | CONTROL RELEASE | 12 | V475G |
| | | 1 | ● | 70 | | 36 | V468H |
| | | 1 | ● | 90 | | 36 | V474H |
| | | 2 | ● | 90 | | 36 | V482H |
|  CP | 40 mm | 2-0 | ● | 70 | | 36 | V478H |
| | | 0 | ● | 70 | | 36 | V479H |
| | | 0 | ● | 90 | | 36 | V485H |
| | | 1 | ● | 70 | | 36 | V480H |
| | | 1 | ● | 90 | | 36 | V486H |
|  CPX | 48 mm | 0 | ● | 90 | | 36 | V1057H |
| | | 1 | ● | 90 | | 36 | V1058H |
| | | 2 | ● | 90 | | 36 | V1059H |

| ▼ 1/2 Circle Reverse Cutting PRIME* Needles | | | | | | | |
|--|-------|-------|-------------|-------------------|----------------|------|---------|
| Needle | Gauge | Color | Length (cm) | Other information | Number per box | Code | |
|  PS-4 | 16 mm | 5-0 | ○ | 45 | Compound curve | 12 | W9436 |
| | | 4-0 | ○ | 45 | MULTIPASS* | 36 | MPV507H |
| | | 4-0 | ○ | 45 | | 12 | V507G |
| | | 4-0 | ○ | 45 | Compound curve | 12 | W9437 |





| ▼ 3/8 Circle Reverse Cutting PRIME* Needles | | | | | | | |
|---|-------|-------|-------------|-------------------|----------------|------|----------|
| Needle | Gauge | Color | Length (cm) | Other information | Number per box | Code | |
|  P-1 | 11 mm | 6-0 | ○ | 45 | MULTIPASS* | 36 | MPV489H |
| | | 6-0 | ○ | 45 | | 36 | V489H |
| | | 5-0 | ○ | 45 | MULTIPASS* | 36 | MPV490H |
| | | 5-0 | ○ | 45 | | 36 | V490H |
|  P-3 | 13 mm | 6-0 | ○ | 45 | MULTIPASS* | 36 | MPV492H |
| | | 6-0 | ○ | 45 | | 36 | V492H |
| | | 5-0 | ○ | 45 | MULTIPASS* | 36 | MPV493H |
| | | 5-0 | ○ | 45 | | 36 | V493H |
| | | 4-0 | ○ | 45 | MULTIPASS* | 36 | MPV494H |
| | | 4-0 | ○ | 45 | | 36 | V494H |
|  PS-3 | 16 mm | 5-0 | ○ | 45 | MULTIPASS* | 36 | MPV500H |
| | | 4-0 | ○ | 45 | | 24 | W9508T |
|  PS-2 | 19 mm | 5-0 | ○ | 45 | MULTIPASS* | 36 | MPV495H |
| | | 4-0 | ○ | 45 | MULTIPASS* | 36 | MPV496H |
| | | 4-0 | ○ | 45 | | 24 | W9515T |
| | | 3-0 | ○ | 45 | MULTIPASS* | 36 | MPV497H |
| | | 3-0 | ○ | 45 | | 36 | V497H |
|  PS-2 S | 19 mm | 4-0 | ○ | 45 | Slimmed | 36 | V496H |
|  PS-1 | 24 mm | 3-0 | ○ | 45 | MULTIPASS* | 36 | MPV683H |
| | | 3-0 | ○ | 45 | | 36 | V683H |
| | | 2-0 | ○ | 45 | MULTIPASS* | 36 | MPV684H |
| | | 2-0 | ○ | 45 | | 36 | V684H |
|  PS | 26 mm | 3-0 | ○ | 70 | MULTIPASS* | 36 | MPV9526H |
|  PSL | 30 mm | 2-0 | ● | 70 | | 36 | V686H |



| ▼ 3/8 Circle Reverse Cutting Needles | | | | | | | |
|--|-------|-------|-------------|-------------------|----------------|------|-------|
| Needle | Gauge | Color | Length (cm) | Other information | Number per box | Code | |
|  C-3 | 13 mm | 6-0 | ● | 45 | | 36 | V384H |
| | | 5-0 | ● | 45 | | 36 | V385H |
| | | 4-0 | ● | 45 | | 36 | V386H |
|  FS-3 | 16 mm | 6-0 | ● | 45 | | 36 | V387H |
| | | 5-0 | ● | 45 | | 36 | V388H |



VICRYL*
(Polyglactin 910) Braided and Monofilament


| ▼ 3/8 Circle Reverse Cutting Needles | | | | | | |
|--|-------|-------|-------------|-------------------|----------------|--------|
| Needle | Gauge | Color | Length (cm) | Other information | Number per box | Code |
|  FS-2 19 mm | 5-0 | ○ | 45 | | 36 | V291H |
| | 5-0 | ● | 45 | | 36 | V391H |
| | 5-0 | ○ | 70 | | 36 | V421H |
| | 4-0 | ○ | 45 | | 36 | V292ZH |
| | 4-0 | ● | 45 | | 12 | V392ZG |
| | 3-0 | ○ | 45 | | 36 | V293H |
| | 3-0 | ● | 45 | | 12 | V393G |
| | 3-0 | ● | 45 | | 36 | V393H |
| | 3-0 | ○ | 70 | | 36 | V423H |
| | 3-0 | ● | 70 | | 36 | V398H |
|  FS-2S 19 mm | 4-0 | ○ | 45 | Slimmed | 36 | V292H |
| | 4-0 | ● | 45 | Slimmed | 36 | V392H |
| | 4-0 | ○ | 70 | Slimmed | 36 | V422H |
| | 4-0 | ● | 70 | Slimmed | 36 | V397H |
|  FS-1 24 mm | 4-0 | ○ | 70 | | 36 | V441H |
| | 4-0 | ● | 70 | | 36 | V451H |
| | 3-0 | ○ | 70 | | 36 | V442H |
| | 3-0 | ● | 70 | | 36 | V452H |
| | 2-0 | ○ | 70 | | 36 | V443H |
| | 2-0 | ● | 70 | | 36 | V453H |
|  FSL 30 mm | 3-0 | ○ | 70 | | 36 | V588H |
| | 3-0 | ● | 70 | | 36 | V585H |
| | 2-0 | ○ | 75 | | 36 | JV589 |
| | 2-0 | ● | 70 | | 36 | V586H |
| | 0 | ○ | 70 | | 36 | V520H |
| | 0 | ● | 70 | | 36 | V587H |

| ▲ 3/8 Circle Conventional Cutting PRIME [®] Needles | | | | | | |
|--|-------|-------|-------------|------------------------|----------------|----------|
| Needle | Gauge | Color | Length (cm) | Other information | Number per box | Code |
|  PC-3 16 mm | 5-0 | ○ | 45 | MULTIPASS [®] | 36 | MPV9505H |
| | 5-0 | ○ | 45 | | 24 | W9505T |
| | 4-0 | ○ | 45 | MULTIPASS [®] | 36 | MPV9506H |
|  PC-5 19 mm | 3-0 | ○ | 70 | MULTIPASS [®] | 36 | MPV9571H |
|  PC-25 26 mm | 4-0 | ○ | 70 | MULTIPASS [®] | 36 | MPV9580H |
| | 4-0 | ○ | 75 | | 24 | W9580T |
| | 3-0 | ○ | 70 | MULTIPASS [®] | 36 | MPV9581H |
| | 3-0 | ○ | 75 | | 24 | W9581T |

| ▲ Straight Cutting Needles | | | | | | | |
|--|-------|-------|-------------|-------------------|----------------|---------|--|
| Needle | Gauge | Color | Length (cm) | Other information | Number per box | Code | |
|  SC-1 13 mm | 4-0 | ● | 45 | | 36 | V383H | |
|  KS-25 25 mm | 3-0 | ● | 70 | | 12 | V2530G | |
|  TS 48 mm | 3-0 | ○ | 70 | | 36 | AB2673H | |
|  KS 60 mm | 4-0 | ○ | 70 | | 36 | V662H | |
| | 3-0 | ● | 70 | | 36 | V523H | |
| | 3-0 | ○ | 70 | | 36 | V663H | |
| | 2-0 | ● | 70 | | 36 | V524H | |
| | 2-0 | ○ | 70 | | 36 | V664H | |




| ▼ 1/2 Circle MICROPOINT* Spatula | | | | | | | |
|---|-------|-------|-------------|-------------------|----------------|-------|--|
| Needle | Gauge | Color | Length (cm) | Other information | Number per box | Code | |
|  GS-14 GS-14 5,5 mm | 10-0 | ● | 30 | | 12 | V450G | |
| | 9-0 | ● | 30 | | 12 | V449G | |
| | 8-0 | ● | 20 | | 12 | V574G | |
|  TG175-8 7 mm | 7-0 | ● | 45 | | 12 | W9565 | |


| ▼ 3/8 Circle MICROPOINT* Spatula | | | | | | | |
|---|-------|-------|-------------|-------------------|----------------|-------|--|
| Needle | Gauge | Color | Length (cm) | Other information | Number per box | Code | |
|  GS-9 6,5 mm | 8-0 | ● | 30 | | 12 | V120G | |
| | 7-0 | ● | 30 | | 12 | V119G | |
|  GS-9 GS-9 6,5 mm | 10-0 | ● | 30 | Monofilament | 12 | V440G | |
| | 9-0 | ● | 30 | | 12 | V549G | |
| | 9-0 | ● | 30 | Monofilament | 12 | V439G | |
| | 8-0 | ● | 20 | | 12 | V547G | |
| | 8-0 | ● | 30 | | 12 | V548G | |
| | 7-0 | ● | 45 | | 12 | V546G | |


| ▼ 1/4 Circle MICROPOINT* Spatula | | | | | | | |
|---|-------|-------|-------------|-------------------|----------------|--------|--|
| Needle | Gauge | Color | Length (cm) | Other information | Number per box | Code | |
|  TG100-8 6,5 mm | 6-0 | ● | 45 | | 12 | JV7544 | |


VICRYL*


(Polyglactin 910) Braided and Monofilament




| ▼ 3/8 Circle MICROPOINT* Reverse Cutting Needles | | | | | | |
|---|-------|-------|-------------|-------------------|----------------|----------|
| Needle | Gauge | Color | Length (cm) | Other information | Number per box | Code |
|  G-6 | 8 mm | 6-0 | ● | 45 | | 12 V551G |
|  G-1 | 11 mm | 6-0 | ● | 45 | | 12 V554G |
|  G-3 | 13 mm | 5-0 | ● | 45 | | 12 V557G |





| ▼ 3/8 Circle CS ULTIMA* Spatula | | | | | | |
|--|--------|-------|-------------|-------------------|----------------|----------|
| Needle | Gauge | Color | Length (cm) | Other information | Number per box | Code |
|  UCS-35 | 6,2 mm | 10-0 | ● | 15 | MONOFILAMENT | 12 V429G |

| ▼ 1/8 Circle CS ULTIMA* Spatula | | | | | | |
|---|--------|-------|-------------|-------------------|----------------|----------|
| Needle | Gauge | Color | Length (cm) | Other information | Number per box | Code |
|  UCS-26 | 6,2 mm | 10-0 | ● | 10 | MONOFILAMENT | 12 V865G |

| ▼ 3/8 Circle Spatula Needles | | | | | | |
|---|-------|-------|-------------|-------------------|----------------|----------|
| Needle | Gauge | Color | Length (cm) | Other information | Number per box | Code |
|  SPAT-10 | 11 mm | 5-0 | ○ | 45 | | 12 W9761 |

| ▼ 1/2 Circle Spatula Needles | | | | | | |
|--|-------|-------|-------------|-------------------|----------------|----------|
| Needle | Gauge | Color | Length (cm) | Other information | Number per box | Code |
|  S-22 | 8 mm | 6-0 | ● | 45 | | 12 W9566 |
| | | 6-0 | ○ | 45 | | 12 W9756 |

| ▼ 1/4 Circle Spatula Needles | | | | | | |
|--|-------|-------|-------------|-------------------|----------------|----------|
| Needle | Gauge | Color | Length (cm) | Other information | Number per box | Code |
|  S-4 | 8 mm | 4-0 | ● | 45 | | 12 V572G |
|  S-29 | 8 mm | 6-0 | ● | 45 | | 12 V555G |
|  S-14 | 8 mm | 6-0 | ● | 45 | | 12 V570G |
| | | 6-0 | ● | 45 | VISI-BLACK* | 12 V575G |
| | | 6-0 | ○ | 45 | | 12 V670G |
| | | 6-0 | ● | 20 | | 12 W9562 |
| | | 5-0 | ● | 45 | VISI-BLACK* | 12 V577G |
| | | 5-0 | ● | 45 | | 12 V571G |
| | | 5-0 | ○ | 45 | | 12 V671G |

| 3/8 Circle Taper Point Needles | | | | | | | |
|--|--------|-------|-------------|-------------------|----------------|------|-------|
| Needle | Gauge | Color | Length (cm) | Other information | Number per box | Code | |
| BV-6  | 4,8 mm | 10-0 | ● | 10 | Monofilament | 12 | V970G |
| BV-4  | 4,9 mm | 10-0 | ● | 15 | Monofilament | 12 | V971G |
| BV130-5  | 6,5 mm | 9-0 | ● | 15 | Monofilament | 12 | V239G |
| | | 8-0 | ● | 20 | | 12 | V400G |
| BV-1  | 9,3 mm | 6-0 | ● | 45 | | 12 | V301G |

| Ligatures | | | | | | | |
|-----------|-------|-------|-------------|-------------------|----------------|--------|--|
| Needle | Gauge | Color | Length (cm) | Other information | Number per box | Code | |
| LIGAPAK* | 4-0 | ● | 250 | Reel | 12 | V1204G | |
| | 3-0 | ● | 250 | Reel | 12 | V1205G | |
| | 2-0 | ● | 250 | Reel | 12 | V1206G | |
| | 0 | ● | 250 | Reel | 12 | V1207G | |
| | 1 | ● | 250 | Reel | 12 | V1208G | |
| | 4-0 | ● | 125 | Reel | 12 | V204G | |
| | 3-0 | ● | 125 | Reel | 12 | V205G | |
| | 2-0 | ● | 125 | Reel | 12 | V206G | |
| | 0 | ● | 125 | Reel | 12 | V207G | |
| SUTUPAK* | 3-0 | ● | 140 | | 36 | V614H | |
| | 2-0 | ● | 140 | | 36 | V615H | |
| | 0 | ● | 140 | | 36 | V616H | |
| | 1 | ● | 140 | | 36 | V617H | |
| | 2 | ● | 140 | | 36 | V618H | |
| SUTUPAK* | 4-0 | ● | 2x70 | | 24 | V623E | |
| | 4-0 | ○ | 2x70 | | 24 | V651E | |
| | 3-0 | ● | 2x70 | | 24 | V624E | |
| | 3-0 | ○ | 2x70 | | 24 | V652E | |
| | 2-0 | ● | 2x70 | | 24 | V625E | |
| | 2-0 | ○ | 2x70 | | 24 | V653E | |
| | 0 | ● | 2x70 | | 24 | V626E | |
| | 0 | ○ | 2x70 | | 24 | V654E | |
| | 1 | ● | 2x70 | | 36 | V627H | |
| | 2 | ● | 2x70 | | 36 | V628H | |

VICRYL*
(Polyglactin 910) Braided and Monofilament

| Ligatures | | | | | | |
|-----------|-------|-------|-------------|-------------------|----------------|--------|
| Needle | Gauge | Color | Length (cm) | Other information | Number per box | Code |
| SUTUPAK* | 5-0 | ● | 3x45 | | 36 | V632H |
| | 5-0 | ○ | 3x45 | | 36 | V642H |
| | 4-0 | ● | 3x45 | | 36 | V633H |
| | 4-0 | ○ | 3x45 | | 36 | V643H |
| | 3-0 | ● | 3x45 | | 36 | V634H |
| | 3-0 | ○ | 3x45 | | 36 | V644H |
| | 2-0 | ● | 3x45 | | 36 | V635H |
| | 2-0 | ○ | 3x45 | | 36 | V645H |
| | 0 | ● | 3x45 | | 24 | V636E |
| | 1 | ● | 3x45 | | 24 | V637E |
| | 1 | ○ | 3x45 | | 24 | V647E |
| | 2 | ● | 3x45 | | 36 | V638H |
| SUTUPAK* | 4-0 | ● | 5x70 | | 24 | V1214E |
| | 3-0 | ● | 5x70 | | 24 | V1215E |
| | 2-0 | ● | 5x70 | | 24 | V1216E |
| | 0 | ● | 5x70 | | 24 | V1217E |
| | 0 | ○ | 5x70 | | 24 | V965E |
| | 1 | ● | 5x70 | | 36 | V1218H |
| | 2 | ● | 5x70 | | 36 | V1219H |
| SUTUPAK* | 5-0 | ● | 6x45 | | 36 | V1223H |
| | 4-0 | ● | 6x45 | | 36 | V1224H |
| | 3-0 | ● | 6x45 | | 36 | V1225H |
| | 3-0 | ○ | 6x45 | | 36 | V978H |
| | 2-0 | ● | 6x45 | | 36 | V1226H |
| | 0 | ● | 6x45 | | 24 | V1227E |
| | 1 | ● | 6x45 | | 24 | V1228E |
| | 2 | ● | 6x45 | | 36 | V1229H |
| SUTUPAK* | 4-0 | ● | 12x45 | | 36 | V903H |
| | 3-0 | ● | 12x45 | | 36 | V904H |
| | 3-0 | ○ | 12x45 | | 36 | V910H |
| | 2-0 | ● | 12x45 | | 24 | V905E |
| | 2-0 | ○ | 12x45 | | 24 | V911E |
| | 0 | ● | 12x45 | | 24 | V906E |
| | 0 | ○ | 12x45 | | 24 | V912E |
| | 1 | ● | 12x45 | | 24 | V907E |