



# SafeStep™

HUBER NEEDLE SET

The Right Step for You and Your Patient

**BAIRD**

## Features and Benefits

### Easy-to-Use:


- ▶ Smooth, simple, single-step safety activation
- ▶ Clear plastic base provides excellent site visibility
- ▶ Needleless access connector for needle-free access
- ▶ Dual port adaptability

### Designed for Patient Comfort:

- ▶ Smallest footprint of all BARD® Safety Infusion Sets
- ▶ Non-absorbant patient comfort pad
- ▶ Improved dressing profile
- ▶ Lower insertion and penetration force than other currently marketed non-siliconized port access needles

### Safety:

- ▶ Built in safety mechanism reduces risk of needle-stick injury



## SAFESTEP™ Huber Needle Set

Codes without "Y" - Injection Site

Product Code	Size
LH0039	19 Ga. x 0.5" (12.7mm)
LH0033	19 Ga. x 0.75" (19mm)
LH0034	19 Ga. x 1" (25mm)
LH0040	19 Ga. x 1.5" (38mm)
LH0037	20 Ga. x 0.5" (12.7mm)
LH0031	20 Ga. x 0.75" (19mm)
LH0032	20 Ga. x 1" (25mm)
LH0038	20 Ga. x 1.5" (38mm)
LH0035	22 Ga. x 0.5" (12.7mm)
LH0029	22 Ga. x 0.75" (19mm)
LH0030	22 Ga. x 1" (25mm)
LH0036	22 Ga. x 1.5" (38mm)

Minimum Gravity Flow Rate for SAFESTEP™ SIS

19 Ga.	1680 ml/hr
20 Ga.	960 ml/hr
22 Ga.	312 ml/hr

Codes with "Y" - Injection Site

Product Code	Size
LH0039YN	19 Ga. x 0.5" (12.7mm)
LH0033YN	19 Ga. x 0.75" (19mm)
LH0034YN	19 Ga. x 1" (25mm)
LH0040YN	19 Ga. x 1.5" (38mm)
LH0037YN	20 Ga. x 0.5" (12.7mm)
LH0031YN	20 Ga. x 0.75" (19mm)
LH0032YN	20 Ga. x 1" (25mm)
LH0038YN	20 Ga. x 1.5" (38mm)
LH0035YN	22 Ga. x 0.5" (12.7mm)
LH0029YN	22 Ga. x 0.75" (19mm)
LH0030YN	22 Ga. x 1" (25mm)
LH0036YN	22 Ga. x 1.5" (38mm)

Priming Volume for SAFESTEP™ SIS

19, 20, and 22 Ga.	with Y-Site, 0.7 ml
19, 20, and 22 Ga.	without Y-Site, 0.3 ml



## Learn More

For more information on Safety Infusion Sets, please contact your Bard Access Systems representative. Or visit us on the web at [www.bardaccess.com](http://www.bardaccess.com)

**Bard France S.A.S.** (France, North Africa) FRANCE  
Tel: + 33 1 39 30 58 58 Fax: + 33 1 39 30 58 59  
service.commande@crbard.com www.bardfrance.fr

**Bard de España S.A.** SPAIN  
Tel: + 34 93 253 7800 Fax: + 34 93 253 7834  
spain.customerservice@crbard.com

**C. R. Bard Portugal L.D.A.** PORTUGAL  
Tel: + 351 21 319 0330 Fax: + 351 21 319 0339  
Bard.Portugal@crbard.com

**Bard S.p.A.** (Italy, Eastern Mediterranean, Middle East) ITALY  
Tel: + 39 06 524 93 249 Fax: + 39 06 529 58 52  
customerservice.italia@crbard.com www.barditalia.it

**Bard Hellas SA** GREECE  
Tel: + 30 210 96 90 780 1 Fax: + 30 210 96 28 810  
Bard.Hellas@crbard.com

**Bard Medica S.A.** SWITZERLAND  
Tel: + 41 44 722 5360 Fax: + 41 44 722 5370  
info-bce.access@crbard.com

**Bard Czech Republic S.R.O.** CZECH REPUBLIC  
Tel: + 420 242 408 630 Fax: + 420 242 410 185  
info-bce.access@crbard.com

**C. R. BARD GmbH** GERMANY  
Tel: + 49 721 9445 124 Fax: + 49 721 9445 100  
info-bce.access@crbard.com

**Bard Poland Sp. z o.o.** POLAND  
Tel: + 48 22 321 0930 Toll Free: 0800 494900  
Fax: + 48 22 321 0938 Toll Free: 0800 494901  
info-bce.access@crbard.com

**Bard Medica S.A.** AUSTRIA  
Tel: + 43 1 49 49 130 Fax: + 43 1 49 49 130 30  
info-bce.access@crbard.com

**Bard Limited** UNITED KINGDOM  
Tel: + 44 1293 529555 Fax: + 44 1293 552428  
customer.services@crbard.com www.barduk.com

**IRELAND**  
Tel: + 353 1 2130706 Fax: + 353 1 2130604  
customer.services@crbard.com

**Bard Norden AB** SWEDEN  
Tel: + 46 42 38 60 00 Fax: + 46 42 38 60 10  
support.nordic@crbard.com

**Bard Norway AS** NORWAY  
Tel: + 47 815 00 266 Fax: + 47 815 00 299  
support.nordic@crbard.com

**Bard Denmark AB** DENMARK  
Tel: + 45 49 26 20 08 Fax: + 45 49 26 20 09  
support.nordic@crbard.com

**Bard Finland Oy** FINLAND  
Tel: + 358 9 347 30 53 Fax: + 358 9 347 10 07  
support.nordic@crbard.com

**Bard Medical SA (Pty) Ltd.** SOUTH AFRICA  
Tel: + 27 11 524 9900 Fax: + 27 86 537 7250  
www.bard.co.za

**Bard Benelux N.V.** BELGIUM  
Tel: + 32 14 28 69 52 Fax: + 32 14 28 69 66  
custservicebelgium@crbard.com

**THE NETHERLANDS**  
Tel: + 31 30 60 00 572 Fax: + 31 30 60 22 345  
custservicenederland@crbard.com