



SafeStep™

HUBER NEEDLE SET

The Right Step for You and Your Patient

BAIRD

Features and Benefits

Easy-to-Use:


- ▶ Smooth, simple, single-step safety activation
- ▶ Clear plastic base provides excellent site visibility
- ▶ Needleless access connector for needle-free access
- ▶ Dual port adaptability

Designed for Patient Comfort:

- ▶ Smallest footprint of all BARD® Safety Infusion Sets
- ▶ Non-absorbant patient comfort pad
- ▶ Improved dressing profile
- ▶ Lower insertion and penetration force than other currently marketed non-siliconized port access needles

Safety:

- ▶ Built in safety mechanism reduces risk of needle-stick injury



SAFESTEP™ Huber Needle Set

Codes without "Y" - Injection Site

| Product Code | Size |
|--------------|------------------------|
| LH0039 | 19 Ga. x 0.5" (12.7mm) |
| LH0033 | 19 Ga. x 0.75" (19mm) |
| LH0034 | 19 Ga. x 1" (25mm) |
| LH0040 | 19 Ga. x 1.5" (38mm) |
| LH0037 | 20 Ga. x 0.5" (12.7mm) |
| LH0031 | 20 Ga. x 0.75" (19mm) |
| LH0032 | 20 Ga. x 1" (25mm) |
| LH0038 | 20 Ga. x 1.5" (38mm) |
| LH0035 | 22 Ga. x 0.5" (12.7mm) |
| LH0029 | 22 Ga. x 0.75" (19mm) |
| LH0030 | 22 Ga. x 1" (25mm) |
| LH0036 | 22 Ga. x 1.5" (38mm) |

Minimum Gravity Flow Rate for SAFESTEP™ SIS

| | |
|--------|------------|
| 19 Ga. | 1680 ml/hr |
| 20 Ga. | 960 ml/hr |
| 22 Ga. | 312 ml/hr |

Codes with "Y" - Injection Site

| Product Code | Size |
|--------------|------------------------|
| LH0039YN | 19 Ga. x 0.5" (12.7mm) |
| LH0033YN | 19 Ga. x 0.75" (19mm) |
| LH0034YN | 19 Ga. x 1" (25mm) |
| LH0040YN | 19 Ga. x 1.5" (38mm) |
| LH0037YN | 20 Ga. x 0.5" (12.7mm) |
| LH0031YN | 20 Ga. x 0.75" (19mm) |
| LH0032YN | 20 Ga. x 1" (25mm) |
| LH0038YN | 20 Ga. x 1.5" (38mm) |
| LH0035YN | 22 Ga. x 0.5" (12.7mm) |
| LH0029YN | 22 Ga. x 0.75" (19mm) |
| LH0030YN | 22 Ga. x 1" (25mm) |
| LH0036YN | 22 Ga. x 1.5" (38mm) |

Priming Volume for SAFESTEP™ SIS

| | |
|--------------------|------------------------|
| 19, 20, and 22 Ga. | with Y-Site, 0.7 ml |
| 19, 20, and 22 Ga. | without Y-Site, 0.3 ml |



Learn More

For more information on Safety Infusion Sets, please contact your Bard Access Systems representative. Or visit us on the web at www.bardaccess.com

Bard France S.A.S. (France, North Africa) FRANCE
Tel: + 33 1 39 30 58 58 Fax: + 33 1 39 30 58 59
service.commande@crbard.com www.bardfrance.fr

Bard de España S.A. SPAIN
Tel: + 34 93 253 7800 Fax: + 34 93 253 7834
spain.customerservice@crbard.com

C. R. Bard Portugal L.D.A. PORTUGAL
Tel: + 351 21 319 0330 Fax: + 351 21 319 0339
Bard.Portugal@crbard.com

Bard S.p.A. (Italy, Eastern Mediterranean, Middle East) ITALY
Tel: + 39 06 524 93 249 Fax: + 39 06 529 58 52
customerservice.italia@crbard.com www.barditalia.it

Bard Hellas SA GREECE
Tel: + 30 210 96 90 780 1 Fax: + 30 210 96 28 810
Bard.Hellas@crbard.com

Bard Medica S.A. SWITZERLAND
Tel: + 41 44 722 5360 Fax: + 41 44 722 5370
info-bce.access@crbard.com

Bard Czech Republic S.R.O. CZECH REPUBLIC
Tel: + 420 242 408 630 Fax: + 420 242 410 185
info-bce.access@crbard.com

C. R. Bard GmbH GERMANY
Tel: + 49 721 9445 124 Fax: + 49 721 9445 100
info-bce.access@crbard.com

Bard Poland Sp. z o.o. POLAND
Tel: + 48 22 321 0930 Toll Free: 0800 494900
Fax: + 48 22 321 0938 Toll Free: 0800 494901
info-bce.access@crbard.com

Bard Medica S.A. AUSTRIA
Tel: + 43 1 49 49 130 Fax: + 43 1 49 49 130 30
info-bce.access@crbard.com

Bard Limited UNITED KINGDOM
Tel: + 44 1293 529555 Fax: + 44 1293 552428
customer.services@crbard.com www.barduk.com

IRELAND
Tel: + 353 1 2130706 Fax: + 353 1 2130604
customer.services@crbard.com

Bard Norden AB SWEDEN
Tel: + 46 42 38 60 00 Fax: + 46 42 38 60 10
support.nordic@crbard.com

Bard Norway AS NORWAY
Tel: + 47 815 00 266 Fax: + 47 815 00 299
support.nordic@crbard.com

Bard Denmark AB DENMARK
Tel: + 45 49 26 20 08 Fax: + 45 49 26 20 09
support.nordic@crbard.com

Bard Finland Oy FINLAND
Tel: + 358 9 347 30 53 Fax: + 358 9 347 10 07
support.nordic@crbard.com

Bard Medical SA (Pty) Ltd. SOUTH AFRICA
Tel: + 27 11 524 9900 Fax: + 27 86 537 7250
www.bard.co.za

Bard Benelux N.V. BELGIUM
Tel: + 32 14 28 69 52 Fax: + 32 14 28 69 66
custservicebelgium@crbard.com

THE NETHERLANDS
Tel: + 31 30 60 00 572 Fax: + 31 30 60 22 345
custservicenederland@crbard.com