
Getinge SMART trolley range
Ergonomic flow of sterile goods

GETINGE GROUP



A close-up photograph of a wheelchair's adjustable armrest and backrest mechanism. The image shows a polished metal armrest curving downwards, attached to a white plastic base. Below the base is a black plastic backrest with a white oval cutout. The background is plain white.

SMART moves

Different people, different positions

The Getinge SMART concept – a whole new level of **convenience**

An ergonomic concept

Loading and distribution systems are an important part of the workflow in operating rooms and sterile supply units. Getinge SMART is a line of carefully designed trolleys that takes loading, unloading and distribution of instruments to a whole new level, both in terms of ergonomics and function.

The trolley range is based on an ergonomic concept for efficient loading and distribution of goods with the highest hygiene and usability. Getinge SMART means top-of-the-line functionality in every detail, from the choice of materials and surfaces to maneuverability, cleanability and ergonomics.

Getinge SMART trolleys can be equipped with a semi-automatic loading function to reduce manual handling of goods when loading or unloading. The trolleys are a joy to handle. The Power-Drive function is an option enabling the user to maneuver the trolley to the desired position with little effort. These optional features assure good ergonomics and give your staff the best possible working conditions.



SMART loading and unloading trolleys

– ergonomic in every detail

Ergonomics is about people being able to work without exertion, strain or injury. It's about moving heavy loads easily. It's about being flexible and understanding that people differ in

height and strength – and in the tasks they perform. Our loading and unloading trolleys* are based on a SMART ergonomic concept that covers every detail of loading and distribution.

SMART fixed height trolleys

- Easy to maneuver
- Easy to clean
- Ergonomic
- SMART loading trolley without “SMART options” can be washed and disinfected in Getinge’s cart washer-disinfectors. Except trolley for HS33.



SMART options

Semi-automatic loading (SAL)

For automatically loading and unloading sterilizer racks, activated by the control integrated in the handle.

SMART adjustable height trolleys (AHT)

Adjust to ergonomic working heights to reduce strain for all users.

Power-Drive (PD)

Significantly reduced strain by enabling moving and maneuvering of trolleys to the desired position with little effort.



*Some variants can be washed and disinfected in Getinge’s cart washer-disinfectors.



A trolley with SMART options has integrated controls in the handle to enable safe use and maneuvering without letting go.

SMART loading and unloading trolleys

technical information

Getinge SMART loading and unloading trolleys* give you the possibility to choose a solution that really meets your needs in every situation.

Loading trolley for rack – Getinge HS33 sterilizers

Sterilizer	Item no	Product description	Max load
HS33	564710370	SMART FHT HS33	30 kg / 66 lbs
HS33+	564710372	SMART FHT HS33+	30 kg / 66 lbs



Loading trolley for rack – Getinge 522/533-HC sterilizers

Sterilizer	Item no	Product description	Max load
522/533-HC	564710471	LT FIX FOR 500	130 kg / 287 lbs



Loading trolley for shelf rack – Getinge 66, 67 and 633HC sterilizers

Full-length load carrier

Sterilizer	Item no	Product description	Max load	Length load carrier
HS/GE 6606, 633-07	6020720670	SMART FHT 66/67x06/13HL/633-07	200 kg / 440 lbs	640 mm / 25.2 in
HS/GE 6617, GSS6717, 633-10	6020720672	SMART FHT 66/67x10/20HL/633-10	200 kg / 440 lbs	990 mm / 38.97 in
HS/GE 6613, GSS6713, 633-13	6020720674	SMART FHT 66/67x13/633-13	200 kg / 440 lbs	1290 mm / 50.78 in
HS/GE 6617, GSS6717	6020720676	SMART FHT 66/67x17	200 kg / 440 lbs	1680 mm / 6.14 in



Half-length load carrier. Two load carriers for full sterilizer load.

Sterilizer	Item no	Product description	Max load	Length load carrier
HS/GE 6613, GSS6713	6020720670	SMART FHT 66/67x06/13HL/633-07	200 kg / 440 lbs	640 mm / 25.2 in
HS/GE 6617, GSS6717	6020720672	SMART FHT 66/67x10/20HL/633-10	200 kg / 440 lbs	840 mm / 33.07 in
HS/GE 6620, GSS6720	6020720672	SMART FHT 66/67x10/20HL/633-10	200 kg / 440 lbs	1290 mm / 50.78 in



With SMART option full-length load carrier

Sterilizer	Item no	Product description	Max load	Length load carrier	AHT	SAL	PD
HS/GE 6610 GSS6710	6020721472	SMART FHT SAL 66/67x10, 633-10	200 kg / 440 lbs	990 mm / 38.97 in		x	
HS/GE 6613 GSS6713	6020721474	SMART FHT SAL 66/67x13, 633-13	200 kg / 440 lbs	1290 mm / 50.78 in		x	
HS/GE 6610 GSS6710	6020723072	SMART AHT 66/67/69x10, 633-10	200 kg / 440 lbs	990 mm / 38.97 in	x		
HS/GE 6613 GSS6713	6020723074	SMART AHT 66/67/69x13, 633-13	200 kg / 440 lbs	1290 mm / 50.78 in	x		
HS/GE 6617 GSS6717	6020723076	SMART AHT 66/67x17	200 kg / 440 lbs	1680 mm / 66.14 in	x		
HS/GE 6610 GSS6710	6020724072	SMART AHT SAL 6x10/633-10	200 kg / 440 lbs	990 mm / 38.97 in	x	x	
HS/GE 6613 GSS6713	6020724074	SMART AHT SAL 6x13/633-13	200 kg / 440 lbs	1290 mm / 50.78 in	x	x	
HS/GE 6610 GSS6710	6020724572	SMART AHT PD 6x10/633-10	200 kg / 440 lbs	990 mm / 38.97 in	x		x
HS/GE 6613 GSS6713	6020724574	SMART AHT PD 6x13/633-13	200 kg / 440 lbs	1290 mm / 50.78 in	x		x
HS/GE 6610 GSS6710	6020725072	SMART AHT SAL PD 6x10/633-10	200 kg / 440 lbs	990 mm / 38.97 in	x	x	x
HS/GE 6613 GSS6713	6020725074	SMART AHT SAL PD 6x13/633-13	200 kg / 440 lbs	1290 mm / 50.78 in	x	x	x

Options (below)

*Containers, shelf racks and basket feeders need to be ordered separately.



With SMART option half-length load carriers. Two load carrier for full sterilizer load.

Options (below)

Sterilizer	Item no	Product description	Max load	Length load carrier	AHT	SAL	PD
HS/GE 6617, GSS6717	6020723072	SMART AHT 66/67/69x10, 633-10	200 kg / 440 lbs	840 mm / 33.07 in	x		
HS/GE 6617, GSS6717	6020724572	SMART AHT PD 6x10/633-10	200 kg / 440 lbs	840 mm / 33.07 in	x		x
HS/GE 6620, GSS6720	6020723072	SMART AHT 66/67/69x10, 633-10	200 kg / 440 lbs	990 mm / 38.97 in	x		
HS/GE 6620, GSS6720	6020724572	SMART AHT PD 6x10/633-10	200 kg / 440 lbs	990 mm / 38.97 in	x		x

With SMART option half-length load carrier. Two load carrier for full sterilizer load.

(Only for docking to automatic loading system, for shorter footprint)

Options (below)

Sterilizer	Item no	Product description	Max load	Length load carrier	AHT	SAL	PD
HS/GE 6617 GSS6717	6020723176	SMART AHT 66/67x17 HL UN-/L	200 kg / 440 lbs	840 mm / 33.07 in	x		

Loading trolley for basket feeder – Getinge 66 and 67 sterilizers

The loading trolley together with a basket feeder enable manual loading of basket into a sterilizer.

Sterilizer	Item no	Product description	Max load	Basket feeder
HS/GE 6606	6020720670	SMART FHT 66/67x06/13HL/633-07	125 kg / 275 lbs*	564613070CE
HS/GE 6610, GSS6710	6020720672	SMART FHT 66/67x10/20HL/633-10	125 kg / 275 lbs*	564613270CE
HS/GE 6613, GSS6713	6020720674	SMART FHT 66/67x13/633-13	125 kg / 275 lbs*	564613470CE
HS/GE 6617, GSS6717	6020720676	SMART FHT 66/67x17	125 kg / 275 lbs*	564613670CE



With SMART option

Options (below)

Sterilizer	Item no	Product description	Max load	Basket feeder	AHT	SAL	PD
HS/GE 6610, GSS6710	6020723072	SMART AHT 66/67/69x10, 633-10	125 kg / 275 lbs*	564613270CE	x		
HS/GE 6610, GSS6710	6020724572	SMART AHT PD 6x10/633-10	125 kg / 275 lbs*	564613270CE	x		x
HS/GE 6613, GSS6713	6020723074	SMART AHT 66/67/69x13, 633-13	125 kg / 275 lbs*	564613470CE	x		
HS/GE 6613, GSS6713	6020724574	SMART AHT PD 6x13/633-13	125 kg / 275 lbs*	564613470CE	x		x
HS/GE 6617, GSS6717	6020723076	SMART AHT 66/67x17	125 kg / 275 lbs*	564613670CE	x		

* The capacity of the basket is limited to 125 kg / 275 lbs.

Loading trolley for automatic basket loader/unloader – Getinge 66 and 67 sterilizers

The loading trolley is used to load/unload automatic loader and unloader for basket.
Not possible to use direct to sterilizer.



With SMART option

Options (below)

Sterilizer with automatic loader/unloader	Item no	Product description	Max load	AHT	SAL	PD
HS/GE 6610, GSS6710	6020727072	SMART AHT BASKET 66/67x10 UN-/L	200 kg / 440 lbs	x		
HS/GE 6613, GSS6713	6020727074	SMART AHT BASKET 66/67x13 UN-/L	200 kg / 440 lbs	x		

Loading trolley for racks – Getinge 69 sterilizers

Load carrier Getinge shelf rack. Possible to dock to sterilizer.

Full-length load carrier

Sterilizer	Item no	Product description	Max load	Length load carrier
HS/GE 6910	6020720872	SMART FHT 6910	200 kg / 440 lbs	990 mm / 38.97 in
HS/GE 6913	6020720874	SMART FHT 6913	200 kg / 440 lbs	1290 mm / 50.78 in
GE 6915	6020720875	SMART FHT 6915	200 kg / 440 lbs	1530 mm / 60.23 in



Half-length load carrier. Two load carrier for full sterilizer load

Sterilizer	Item no	Product description	Max load	Length load carrier
HS/GE 6917	6020720872	SMART FHT 6910	200 kg / 440 lbs	840 mm / 33.07 in



With SMART option full-length load carrier

Options (below)

Sterilizer	Item no	Product description	Max load	Length load carrier	AHT	SAL	PD
HS/GE 6910	6020721672	SMART FHT 6910	200 kg / 440 lbs	990 mm / 38.97 in		x	
HS/GE 6913	6020721674	SMART FHT 6913	200 kg / 440 lbs	1290 mm / 50.78 in		x	
HS/GE 6910	6020723072	SMART AHT 66/67/69x10, 633-10	200 kg / 440 lbs	990 mm / 38.97 in	x		
HS/GE 6913	6020723074	SMART AHT 66/67/69x13, 633-13	200 kg / 440 lbs	1290 mm / 50.78 in	x		
GE 6915	6020723075	SMART AHT 6915	210 kg / 462 lbs	1530 mm / 60.23 in	x		
HS/GE 6910	6020724072	SMART AHT SAL 6x10/633-10	200 kg / 440 lbs	990 mm / 38.97 in	x	x	
HS/GE 6913	6020724074	SMART AHT SAL 6x13/633-13	200 kg / 440 lbs	1290 mm / 50.78 in	x	x	
HS/GE 6610, GSS6710	6020724572	SMART AHT PD 6x10/633-10	200 kg / 440 lbs	990 mm / 38.97 in	x		x
HS/GE 6613, GSS6713	6020724574	SMART AHT PD 6x13/633-13	200 kg / 440 lbs	1290 mm / 50.78 in	x		x
GE 6915	6020724575	SMART AHT PD 6915	210 kg / 462 lbs	1530 mm / 60.23 in	x		x
HS/GE 6910	6020725072	SMART AHT SAL PD 6x10/633-10	200 kg / 440 lbs	990 mm / 38.97 in	x	x	x
HS/GE 6913	6020725074	SMART AHT SAL PD 6x13/633-13	200 kg / 440 lbs	1290 mm / 50.78 in	x	x	x
GE 6915	6020725075	SMART AHT SAL PD 6915	210 kg / 462 lbs	1530 mm / 60.23 in	x	x	x

With SMART option half-length load carrier. Two load carrier for full sterilizer load.

Sterilizer	Item no	Product description	Max load	Length load carrier	AHT	SAL	PD
HS/GE 6917	6020723072	SMART AHT 66/67/69x10, 633-10	200 kg / 440 lbs	840 mm / 33.07 in	x		



Load carrier Getinge shelf rack, loading height 1100 mm / 40.30 in. Possible to dock to sterilizer.

With SMART Option Full length Load carrier

Sterilizer	Item no	Product description	Max load	Length load carrier	AHT	SAL	PD
GEB 6910	6020726072	SMART AHT GEB 6910	200 kg / 440 lbs	990 mm / 38.97 in	x		
GEB 6913	6020726074	SMART AHT GEB 6913	200 kg / 440 lbs	1290 mm / 50.78 in	x		
GEB 6915	6020726075	SMART AHT GEB 6915	210 kg / 462 lbs	990 mm / 38.97 in	x		

Loading trolleys for washer-disinfectors

Load carrier Getinge wash cart*. Possible to dock to 8666/8668/88T/CM320 washer-disinfectors and loading equipment for Automation System 2.0.

Single-end loading

Washer-disinfector	Item no	Product description	Max load
8666/8668/88T/CM320	6020720171	SMART FHT WD 8666/8668/88T/CM320	165 kg / 363 lbs

Double-end loading

Washer-disinfector	Item no	Product description	Max load
8666/8668/88T/CM320	6020720172	SMART FHT WD DH 8666/8668/88T/CM320	165 kg / 363 lbs



With SMART option single-end loading

Washer-disinfector	Item no	Product description	Max load	AHT	SAL	PD
8666/8668/88T/CM320	6020720271	SMART AHT WD 8666/8668/88T/CM320	165 kg / 363 lbs	x		

Options (below)

*Wash cart needs to be ordered separately.



SMART distribution trolleys **for ergonomic transport efficiency**

Getinge SMART distribution trolleys are designed for the effective, ergonomic transport of baskets and/or containers. They are used in

packing areas in sterile zones, and for transport to and from soiled zones throughout the hospital.



Distribution trolley hard-cover

In this trolley, instrument containers and baskets can be grouped together with utensils, linen, disposables and anesthesia material for transport to an internal OR department or to an external hospital or storage facility. The hard cover protects the sterile goods from dust and damage during transport.

- Simple and safe transport
- Ergonomic
- Can be washed and disinfected in Getinge cart washer-disinfector.



Distribution trolley open

This trolley is a mobile storage unit for modular baskets and containers. The open design makes viewing and picking of individual packs very easy. Shelves and baskets are easily secured for safe transport.

- Easy to load and unload
- Saves storage space
- Silent and easy to maneuver
- Can be washed and disinfected in Getinge cart washer-disinfector.



The doors can be opened 260° and fixed in the open position – very convenient when loading and unloading. The door handle can be fitted with an optional lock. Shelf supports can be moved, which maximizes utilization of the trolley.



A maximum of five shelves can be used. Three are supplied as standard. Additional shelves are available (optional) for extra storage space. Each shelf is supplied with two dividers to keep goods in place during transport.



To facilitate maneuvering, the trolley is equipped with two lockable swivel casters and two rigid casters. The caster diameter is 160 mm. An optional pin and tow bar is available to enable towing.



The flexible placement of support guides and shelves maximizes trolleys utilization.



Trolleys can be loaded and unloaded from both sides. Shelves and baskets can be pulled out halfway for easier access.



The high-quality swivel casters (two of which are lockable) enable easy to moving and maneuvering even when the trolley is fully loaded.

SMART distribution trolleys

technical information

With Getinge SMART distribution trolleys you always get safe, efficient distribution. The hard-cover model protects sterile goods from dust and damage. The open trolley model gives you a mobile unit for modular baskets and containers.



Distribution trolley hard-cover (DTHC)

Item no	Capacity	Options (Item no)	Max load	Trolley dimensions (LxWxH)
564719070	6 or 8 SPRI/ISO high* or 15 DIN StU ½ size**	564719071 wire shelf incl. 2 dividers 564719075 pin and towbar	145 kg / 319 lbs	1286x699x1502 mm / 50.63x27.52x59.13 in

*SPRI high: 585 × 395 × 195 mm / 23.03 × 15.55 × 7 × 67 in. ISO high: 600 × 400 × 195 mm / 23.62 × 15.75 × 7 × 67 in (for 8 SPRI/ISO, order one extra shelf)

** DIN StU ½ size: 575 × 280 × 135 mm / 22.64 × 11.02 × 5.31 in



DTO double (shelves)
Item No: 564718970



DTO single (mixed)
Item No: 564718370



DTO single SPRI (5 high)
Item No: 564718170
(Baskets not included)



DTO single (shelves)
Item No: 564718570
(Container not included)

Getinge SMART open distribution trolleys are available in single or double versions of each model:

- Three models configured for standard loads. Delivered with guides for baskets, shelves for containers or mixed
- One base model that can be individually configured. For the base model, guides are specified and ordered separately for customers' individual configuration of a mixed load of shelves and/or baskets

Distribution trolley open (DTO)

DTO Single

Type of load	Description	Item no	Max load	Dimensions (L×W×H)
Baskets	DTO single SPRI, 10 low or 5 high	564718170	85 kg / 187 lbs	736 × 429 × 1577 mm / 28.98 × 16.89 × 62.09 in
Mixed	DTO single, 2 SPRI high/3 shelves	564718370	85 kg / 187 lbs	736 × 429 × 1577 mm / 28.98 × 16.89 × 62.09 in
Shelves	DTO single, 5 shelves	564718570	85 kg / 187 lbs	736 × 429 × 1577 mm / 28.98 × 16.89 × 62.09 in
Individually configured	DTO single base without guides	564718070	85 kg / 187 lbs	736 × 429 × 1577 mm / 28.98 × 16.89 × 62.09 in

DTO Double

Type of load	Description	Item no	Max load	Dimensions (L×W×H)
Baskets	DTO single SPRI, 20 low or 10 high	564718770	170 kg / 374 lbs	1416 × 429 × 1577 mm / 55.75 × 16.89 × 62.09 in
Mixed	DTO single, 4 SPRI high/6 shelves	564718870	170 kg / 374 lbs	1416 × 429 × 1577 mm / 55.75 × 16.89 × 62.09 in
Shelves	DTO single, 10 shelves	564718970	170 kg / 374 lbs	1416 × 429 × 1577 mm / 55.75 × 16.89 × 62.09 in
Individually configured	DTO double base without guides	564718670	170 kg / 374 lbs	1416 × 429 × 1577 mm / 55.75 × 16.89 × 62.09 in

Options for Single and Double

	Description	Item no	Max load	Dimensions (L×W×H)
	Pair of guides SPRI	564718171	18 kg / 39 lbs	30 × 425 × 50 mm / 1.18 × 16.73 × 1.99 in
	Pair of guides ISO	564718771	18 kg / 39 lbs	30 × 425 × 50 mm / 1.18 × 16.73 × 1.99 in
	Pair of guides and shelf	564718571	18 kg / 39 lbs	30 × 425 × 50 mm / 1.18 × 16.73 × 1.99 in
	External transport lock	564718470		
	Rack for disposable covers	560181301		
	Disposable transparent covers, 50 pcs/roll	560125070		
	Reusable cover for DTO single	6003000052		
	Reusable cover for DTO double	6003000053		

GETINGE GROUP

Getinge Infection Control AB
P O Box 69
SE-305 05 Getinge, Sweden

www.getingegroup.com

MAQUET
GETINGE GROUP

ARJOHUNTLEIGH
GETINGE GROUP

GETINGE
GETINGE GROUP

Getinge Group is a leading global provider of innovative solutions for operating rooms, intensive-care units, hospital wards, sterilization departments, elderly care and for life science companies and institutions. With a genuine passion for life we build quality and safety into every system. Our unique value proposition mirrors the continuum of care, enhancing efficiency throughout the clinical pathway. Based on our first-hand experience and close partnerships, we are able to exceed expectations from customers – improving the every-day life for people, today and tomorrow.