
GE Healthcare TruSignal sensors and cables – monitor compatibility



TruSignal Integrated Ear sensor with Datex connector TS-E4-N 4 m/13 ft

*Trademark of General Electric Company

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Representative for the most current information.



TruSignal sensors and cables – patient connections

(For more information refer to the sensor's instructions for use)

Patient weight ranges and application suggestions

Preferably apply the sensor to a site where the light source and detector are aligned

TruSignal Adult/Pediatric sensor						20-50 kg (44-110 lb)	≥ 50 kg (≥ 110 lb)
TS-AP-10 TS-AP-25	10/box 25/box						
TruSignal Adult Adhesive Wrap sen	sor					>20	0 kg (>44 lb)
TS-AAW-10 TS-AAW-25	10/box 25/box						
TruSignal Pediatric Adhesive Wrap sensor					3-20 kg (6.6-44 lb)		
(^{**} ^{**} ^{**} ^{**} ^{**} ^{**} ^{**} ^{**}	10/box 25/box						
TruSignal AllFit Sensor		Attachment Accessori	es	≤ 3 kg (≤ 6.6 lb)	3-20 kg (6.6-44 lb)	20-50 kg (44-110 lb)	≥ 50 kg (≥ 110 lb)
TS-AF-10 TS-AF-25 Replacement	10/box 25/box tape included	** : • ** : •	OXY-RT 100/box Replacement adhesive tape OXY-RTB 100/box Replacement adhesive tape - teddy bears	(TruTrak/TruTrak+* and TruSignal monitors only)			
TruSignal Sensitive Skin Sensor		Attachment Accessori	es	≤ 3 kg (≤ 6.6 lb)	3-20 kg (6.6-44 lb)	20-50 kg (44-110 lb)	≥ 50 kg (≥ 110 lb)
	3/box nd adhesive tape included	OXY-SND 3/box Infant foam sandal	OXY-RWS, OXY-RWM 24/box Replacement small and medium foam wrap CXY-RTW 100/box Replacement wide adhesive tape	(TruTrak/TruTrak+ and TruSignal monitors only) \mathcal{M}			
TruSignal Wrap Sensor		Attachment Accessorie	es		3-20 kg (6.6-44 lb)	20-50 kg (44-110 lb)	≥ 50 kg (≥ 110 lb)
	14-H 1/box nd adhesive tape included		OXY-RWL 24/box Replacement large foam wrap OXY-RTW 100/box Replacement wide adhesive tape				
TruSignal Finger Sensor						2	20 kg (≥ 44 lb)
	F-D, TS-F2-GE, TS-F4-GE, F4-MC, TS-F1-H, F4-H, TS-F4-N 1/box						
TruSignal FingerTip Sensor						>	30 kg (> 66 lb)
	SA-D, TS-SA4-GE, SA3-MC 1/box						
TruSignal PediTip Sensor							5-30 kg (33-66 lb)
TS-SP-D, TS-SP3-GE, TS-SP3-MC 1/box				M			
TruSignal Ear Sensor		Attachment Accessori	es		≥ 10 kg (≥ 22 lb)		
TS-E4-N, TS-E	-GE, TS-E4-GE, TS-E4-MC, 4-H 1/box d bedsheet clip included		OXY-HB 3/box Replacement Headband	(Z		

Sensor selection considerations

Older pediatric through adult application

- Medium- to long-term monitoring
- Recommended for:
- high-motion environments
- infection control
- low-perfusion conditions

Older pediatric through adult application

- Medium- to long-term monitoring
- Recommended for:
- high-motion environments
- infection control
- low-perfusion conditions

Infant through older pediatric application

- Medium- to long-term monitoring
- Recommended for:
- high-motion environments
- infection control
- low-perfusion conditions
- patients with fragile skin

Infant through adult application

- Medium- to long-term monitoring
- Recommended for:
- high-motion environments
- infection control
- patients with long fingernails
- arthritic patients

Premature infant/neonate through small adult application

- Medium- to long-term monitoring
- Recommended for:
- high-motion environments if used with adhesive tape
- fragile skin if used with foam wrap or infant foam sandal
- patients with long fingernails

Infant/neonate through adult application

- Short- to medium-term monitoring
- Recommended for:
- high-motion environments if used with adhesive tape
- fragile skin if used with foam wrap or infant foam sandal
- patients with long fingernails

Older pediatric through adult application

- Short-term monitoring
- Recommended for:
- spot-checks

Older pediatric through adult application

• Short- to long-term monitoring

Pediatric through small adult application

• Short- to long-term monitoring

Pediatric through adult application

- Short- to medium-term monitoring
- Recommended for:
- use during stress testing with use of headband
- low-perfusion conditions
- poor peripheral circulation