

699600


Universal Set Standard

Product details

Descriptive feature: Not Applicable

Sterility: Sterile

Images



Products included in set

Product ID	Quantity	Product name	Size	Descriptive feature
750105	1,0	Surgical Gown CLASSIC	L-L	Antistatic treated / Alcohol repellent, Standard performance, Two towels, Wrapped
996106	1,0	Mayo Stand Cover	79cm x 145cm	Absorbent
26153-20	4,0	Cellulose Towel	18cm x 25cm	
910200	2,0	Adhesive OP-Towel	90cm x 75cm	3-ply, Adhesive, Patched
750115	1,0	Surgical Gown CLASSIC	XL-L	Antistatic treated / Alcohol repellent, Standard performance, Two towels, Wrapped

Product ID	Quantity	Product name	Size	Descriptive feature
911500	1,0	OP-Tape	9cm x 49cm	Adhesive
967772	1,0	Adhesive OP-Sheet	175cm x 175cm	3-ply, Adhesive, Patched
967773	1,0	Adhesive OP-Sheet	240cm x 150cm	3-ply, Adhesive, Patched
90190	1,0	Table Cover	150cm x 190cm	Absorbent

Delivered items

699600-07

Country of origin: Belgium

Shelf life: 5 years

Sterilization method: Beta

Packing information: First packaging layer is a peel-open sterile barrier, plastic/plastic. Once opened the sterile barrier cannot be closed again. Second layer is a corrugated board dispenser box. Third layer is a corrugated board transport box.

Packaging level	Pack count	GS1 code
Consumer pack	1	7332430984234
Dispenser box	4	7332430984227
Transport box	8	7332430984210
Pallet	96	7332430984203

Material

Animal tissues:

No

Human blood derivatives:

No

Medicinal substances:

No

Phthalates:

No

Polyvinyl chloride:

No

Technical

Disposal instructions

Non-hazardous waste used BARRIER products and sterility barriers should, in the majority of cases, be classified as non-hazardous waste. They contain high amounts of energy and are well suited for incineration. BARRIER products do not contain any hazardous substances that can leach out if the products are land filled. Transport boxes are designed to fit existing recovery systems. The new BARRIER packaging system complies with the Packaging Waste Directive of the European Union.

Storage instructions

Mölnlycke Health Care recommends that BARRIER products are stored under normal storage conditions. All layers of packaging should be kept intact until access to the underlying layers is needed. Storage facilities for products only protected by the sterility barrier should be kept under conditions where low level of particulate air contamination prevail, so that it would not constitute a risk to the patient when the package is opened and the product is used.

Classification

Regulation type	MDD Class IS
CE Certificate Number:	CE 01966
Notified body medical devices/PPE:	BSI (0086)
Intended use MDD:	Surgical drapes, when sterilised, are intended to minimize the spread of micro-organisms, in order to reduce the risk for post operative wound infection.

Applied standards: The standards presented below is a selection of the most essential standards that are adhered to.

EN 1041, EN 556-1, EN 13795, EN 62366, EN ISO 9001, EN ISO 13485, EN ISO 10993-1, EN ISO 11607-1, EN ISO 11607-2, EN ISO 15223-1, ISO 15223-2, ISO 14001

Removable Label

Yes

GMDN Code (Global Medical Device Nomenclature)

33961