

Ear Tampons


Ear tampons are used for post otoplasty tamponing of the outer acoustic canal and outer ear drum. They are made of diverse materials, most commonly used today are; gelatin, PVA (Polyvinyl alcohol) and cotton.

Ear tampons protect the operated area from external, mechanical influences and can be used in the application of antibiotics.

Tampons, PVA, sterile

- Expanding ear tampon for post-tympanoplastic tamponing of the external acoustic canal
- With or without central lumen
- Single, sterile packaging

ATTENTION:
Not suitable for application in the middle ear or implantation!
Must be kept constantly moist for the entire duration of use!



Art. No.	Size	Feature	Packaging
66400	9 x 15 mm	with lumen	50 pieces / box
66410	12 x 24 mm	without lumen	20 pieces / box
66420	9 x 15 mm	without lumen	50 pieces / box



In the picture above: delivery condition (dry)

In the picture below: expanded state (wet)