

Brand:



<b>Description:</b>	O.R. Towel
<b>Final dimensions:</b>	45 cm x60 cm
<b>Composition:</b>	Polyamide, Polyester, Viscose/Polyester Nonwoven.
<b>Construction:</b>	Texart T
<b>Latex:</b>	Free
<b>PVC:</b>	Free
<b>Re use:</b>	No
<b>Counting Efficiency approx</b>	99,5 % (see note below)
<b>Colour:</b>	White
<b>Sterilized by:</b>	Moist heat
<b>Shelf life:</b>	60 months
<b>CFU/g</b>	N.A.
<b>Size variation</b>	±3 cm
<b>Liquid absorbency time</b>	<2 s
<b>Water soluble substances</b>	<0,5 %
<b>Fluorescence</b>	conform
<b>pH (acidity or alkalinity)</b>	7±1
<b>Ether soluble substances</b>	<0,5 %
<b>Surface-active substances</b>	conform
<b>Extract colouring matter</b>	conform
<b>Dry Linting (&gt;0,5µm)</b>	<5000 particles
<b>Absorption EN 1644-1</b>	780%
<b>Absorption EN 1644-2</b>	10 g/g
<b>Absorption EN 13726-1</b>	13 g/cm <sup>2</sup>
<b>Dry Tensile strength MD</b>	20 N/5cm
<b>Dry Tensile strength TD</b>	120 N/5 cm
<b>Dry bursting strength</b>	250 kPa
<b>Weigth/sqm</b>	120±10% g/m <sup>2</sup>



<b>Nr and Revision Date:</b>	2 -- 2017/03/02
<b>CE Class:</b>	I Steril
<b>CE mark:</b>	CE0120

<b>Directive:</b>	MDD 93/42/ EEC amended by 2007/47/EC
<b>Applicable norms:</b>	EN ISO 13485:2012 ; EN ISO 9001:2008 ; EN 980:2008 ; ISO 15223-1:2012 ; EN 1041+A1:2013 ; EN ISO 11607-1:2009 ; EN ISO 11607-2:2006 ; EN 556-1:2001 ; EN 1644-1:1997 ; EN 1644-2:2000 ; ISO 9073-1:1989 ; ISO 9073-3:1989 ; ISO 9073-6:2000 ; ISO 9073-10:2003 ; EN ISO 17665-1 :2006 ; ISO/TS 17665-2:2009 ; EN ISO 14971:2012 ; EN ISO 10993-1:2009 ; EN ISO 10993-1/AC:2010 ; EN ISO 10993-5:2009 ; EN ISO 10993-18:2009 ; EN 556-1:2001/AC:2006 ;

Level	Type	Dimension	Stamp	Labels
1	2 Peel-pack	160x290mm		labeled
2	2 Peel-pack	200x340mm	stamped	
3	80 Folding box	375x285x450mm		labeled
4	160 Carton box	590x390x465mm		labeled

<b>Important Notes:</b>	Crepe texture with X-ray contrast thread inside Note: Above parameters mean target values
-------------------------	--

<b>Counting Efficiency approx</b>	Surgeons and nurses shall never assume that number of surgical sponges, in each pack, is exact at all times.  The usual counting efficiency of Bastos Viegas products in double packing with control TAG or tied 5's is 99.5%.  The tracking of sponges is full responsibility of the healthcare professionals. In order to avoid retained surgical sponges following surgical operation, (1) follow rigorous manual counting protocols,
-----------------------------------	---

(2) use sponges marked with a radiopaque marker, as well as, (3) careful examination of the abdomen by the operating surgeon before closure.

<b>Storage instructions</b>	Store in a dry place at room temperature and protect from sunlight. The product has good stability and therefore, normal conditions of medical devices storage shall apply
<b>Waste instructions</b>	No special requirements. The waste must be disposed according to local legal regulations for the disposable of medical waste.

<b>WASTE AND ENVIRONMENT REGULATIONS</b>	Weight (gr) of packing material used in 100 pcs	1318
--	---	------

**ISSUED BY: Quality Management Department, 2017/09/13**

---

CONFIDENTIAL: The information contained herein is confidential and belongs to Bastos Viegas, s.a. It is provided for specific purposes and in connection with established agreements. It is not permitted disclosure, reproduction, copying or distribution, in whole or in part, for other purposes or to any third parties.

---