

Brand:



Description:	Abdominal swab
Final dimensions:	30 cm x45 cm
Composition:	Polyamide, Polyester, Viscose/Polyester Nonwoven.
Construction:	Texart C
Latex:	Free
PVC:	Free
Re use:	No
Counting Efficiency approx	99,5 % (see note below)
Colour:	White
Sterilized by:	Ethylene Oxide
Shelf life:	60 months
CFU/g	N.A.
Size variation	±3 cm
Liquid absorbency time	<2 s
Water soluble substances	<0,5 %
Fluorescence	conform
pH (acidity or alkalinity)	7±1
Ether soluble substances	<0,5 %
Surface-active substances	conform
Extract colouring matter	conform
Dry Linting (>0,5µm)	<5000 particles
Weigth/sqm	105±10% g/m ²
Absorption EN 1644-1	1000%
Absorption EN 1644-2	16,5 g/g
Absorption EN 13726-1	18 g/cm ²
Dry Tensile strength MD	120 N/5cm
Dry Tensile strength TD	30 N/5 cm
Dry bursting strength	250 kPa

STERILE EO



isento
Latex
free

isento
PVC
free



Nr and Revision Date:	
CE Class:	Ila
CE mark:	CE0120
Directive:	MDD 93/42/ EEC amended by 2007/47/EC
Applicable norms:	EN ISO 13485:2012 ; EN ISO 9001:2008 ; EN 980:2008 ; ISO 15223-1:2012 ; EN 1041+A1:2013 ; EN ISO 11607-1:2009 ; EN ISO 11607-2:2006 ; EN 556-1:2001 ; EN 1644-1:1997 ; EN 1644-2:2000 ; ISO 9073-1:1989 ; ISO 9073-3:1989 ; ISO 9073-6:2000 ; ISO 9073-10:2003 ; EN ISO 17665-1 :2006 ; ISO/TS 17665-2:2009 ; EN ISO 14971:2012 ; EN ISO 10993-1:2009 ; EN ISO 10993-1/AC:2010 ; EN ISO 10993-5:2009 ; EN ISO 10993-18:2009 ; EN 556-1:2001/AC:2006 ;

Level	Type	Dimension	Stamp	Labels
1	5 Peel-pack	160x235mm		labeled
2	5 Peel-pack	200x290mm	stamped	
3	180 Folding box	375x285x450mm		labeled
4	360 Carton box	590x390x465mm		labeled

Important Notes:	Crepe texture with X-ray contrast thread inside After sterilization there is a shrinking effect and basis weight is approx. 120±10% g/m ² Note: Above parameters mean target values
-------------------------	--

Counting Efficiency approx	Surgeons and nurses shall never assume that number of surgical sponges, in each pack, is exact at all times. The usual counting efficiency of Bastos Viegas products in double packing with control TAG or tied 5's is 99.5%. The tracking of sponges is full responsibility of the healthcare professionals. In order to avoid retained surgical sponges following surgical operation, (1) follow rigorous manual counting protocols,
-----------------------------------	---

(2) use sponges marked with a radiopaque marker, as well as, (3) careful examination of the abdomen by the operating surgeon before closure.

Storage instructions	Store in a dry place at room temperature and protect from sunlight. The product has good stability and therefore, normal conditions of medical devices storage shall apply
Waste instructions	No special requirements. The waste must be disposed according to local legal regulations for the disposable of medical waste.

WASTE AND ENVIRONMENT REGULATIONS	Weight (gr) of packing material used in 100 pcs	499
--	---	-----

ISSUED BY: Quality Management Department, 2017/09/13

CONFIDENTIAL: The information contained herein is confidential and belongs to Bastos Viegas, s.a. It is provided for specific purposes and in connection with established agreements. It is not permitted disclosure, reproduction, copying or distribution, in whole or in part, for other purposes or to any third parties.
