



Product Name Surgifix Zerolatex

Producer
FRA Production S.p.A.

Expiry Date
6 years

Product Description
Elastic netting made of polyamide and elastomeric yarns for medical use.
Available in different sizes.
Class I medical device - CND code: M03030103.

Intended use
Elastic tubular netting for fixing compresses, gauzes and medications.

Reference Regulations
European Directives 93/42/EC and 2007/47/EC and successive amendments.

Certification
UNI EN ISO 9001:2015. ISO 13485:2016. CE MARKING

Warnings
Do not apply directly in contact with the wound. This product is not sterile.
Do not put in contact with flames. Do not swallow. Avoid improper use. Keep out of reach of children.
This product is to be worn, not wrapped.
This product does not contain latex.
If the product is too tight, use a larger size.

Instructions for use

Wear this product in order to keep sanitary medications in place, choosing the appropriate size.

Product storage

Keep away from heat, humidity and light sources.

Related enclosed documents

Declaration of conformity to the reference regulations of each batch. Bill of lading.

Packaging information

Available in different lengths, as specified in the list of the available product versions.

Disposal information

Disposal must be carried out in compliance with related regulations.
In Italy, this product is disposable as "waste from processed textile fibers".

Guaranteed performances

Dimensional characteristics of the product, when stretched, are guaranteed.
This product lets skin breathe and allows maximum freedom of movement, though it keeps the dressing in place.
Biocompatibility tests have been carried out on the product, in order to prove its use suitability.

Composition

86% POLIAMMIDE 14% ELASTOMERO



PRODUCTION S.p.A.

SP000P21
"Surgifix Zerolatex" elastic netting

Revision 011
Date 19/11/19

Product available versions																
Item code	Item description											Sales pack			Number of packs per export cartons	
00250001	SURGIFIX ZEROLATEX TIPO A 25 MT											Mt 25,00			68	

Product technical data																
Weight when stretched (g/m) (due to an applied load of of 200g per counter stitch)			Width (mm) when unstretched			Expansion (mm) due to an applied load of 13Kg			Width (mm) after the expansion			Longitudinal elongation (%) (due to a 200g applied load)			Number of counter stitches	
std	min	max	std	min	max	std	min	max	std	min	max	std	min	max		
0,83	0,74	0,91	13,0	11,7	14,3	44,0	40,0	48,0	12,0	10,8	13,2	290,00	275,00	304,00	022	

Product available versions																
Item code	Item description											Sales pack			Number of packs per export cartons	
00250002	SURGIFIX ZEROLATEX TIPO B 25 MT											Mt 25,00			68	

Product technical data																
Weight when stretched (g/m) (due to an applied load of of 200g per counter stitch)			Width (mm) when unstretched			Expansion (mm) due to an applied load of 13Kg			Width (mm) after the expansion			Longitudinal elongation (%) (due to a 200g applied load)			Number of counter stitches	
std	min	max	std	min	max	std	min	max	std	min	max	std	min	max		
1,57	1,41	1,73	25,0	22,5	27,5	146,0	132,0	160,0	30,0	27,0	33,0	290,00	275,00	304,00	046	

Product available versions																
Item code	Item description											Sales pack			Number of packs per export cartons	
00250003	SURGIFIX ZEROLATEX TIPO C 25 MT											Mt 25,00			52	

Product technical data																
Weight when stretched (g/m) (due to an applied load of of 200g per counter stitch)			Width (mm) when unstretched			Expansion (mm) due to an applied load of 13Kg			Width (mm) after the expansion			Longitudinal elongation (%) (due to a 200g applied load)			Number of counter stitches	
std	min	max	std	min	max	std	min	max	std	min	max	std	min	max		
2,32	2,08	2,55	38,0	34,2	41,8	246,0	221,4	270,6	48,0	43,2	52,8	290,00	275,00	304,00	070	

Product available versions																
Item code	Item description											Sales pack			Number of packs per export cartons	
00250004	SURGIFIX ZEROLATEX TIPO D 25 MT											Mt 25,00			44	

Product technical data																
Weight when stretched (g/m) (due to an applied load of of 200g per counter stitch)			Width (mm) when unstretched			Expansion (mm) due to an applied load of 13Kg			Width (mm) after the expansion			Longitudinal elongation (%) (due to a 200g applied load)			Number of counter stitches	
std	min	max	std	min	max	std	min	max	std	min	max	std	min	max		
3,04	2,74	3,35	49,0	44,1	53,9	350,0	315,0	385,0	76,0	69,0	83,0	290,00	275,00	304,00	094	

Product available versions																
Item code	Item description											Sales pack			Number of packs per export cartons	
00250006	SURGIFIX ZEROLATEX TIPO F 25 MT											Mt 25,00			32	

Product technical data																
Weight when stretched (g/m) (due to an applied load of of 200g per counter stitch)			Width (mm) when unstretched			Expansion (mm) due to an applied load of 13Kg			Width (mm) after the expansion			Longitudinal elongation (%) (due to a 200g applied load)			Number of counter stitches	
std	min	max	std	min	max	std	min	max	std	min	max	std	min	max		
4,63	4,17	5,09	78,0	70,2	85,8	560,0	504,0	616,0	140,0	126,0	154,0	290,00	275,00	304,00	146	

Product available versions																
Item code	Item description											Sales pack			Number of packs per export cartons	
00250007	SURGIFIX ZEROLATEX TIPO G 25 MT											Mt 25,00			24	

Product technical data																
Weight when stretched (g/m) (due to an applied load of of 200g per counter stitch)			Width (mm) when unstretched			Expansion (mm) due to an applied load of 13Kg			Width (mm) after the expansion			Longitudinal elongation (%) (due to a 200g applied load)			Number of counter stitches	
std	min	max	std	min	max	std	min	max	std	min	max	std	min	max		
5,41	4,87	5,96	88,0	80,0	96,0	664,0	597,0	730,0	200,0	180,0	220,0	290,00	275,00	304,00	170	

Product available versions																
Item code	Item description											Sales pack			Number of packs per export cartons	
00250008	SURGIFIX ZEROLATEX TIPO H 25 MT											Mt 25,00			20	

Product technical data																
Weight when stretched (g/m) (due to an applied load of of 200g per counter stitch)			Width (mm) when unstretched			Expansion (mm) due to an applied load of 13Kg			Width (mm) after the expansion			Longitudinal elongation (%) (due to a 200g applied load)			Number of counter stitches	
std	min	max	std	min	max	std	min	max	std	min	max	std	min	max		
7,00	6,30	7,70	113,0	101,2	124,3	850,0	765,0	935,0	250,0	225,0	275,0	290,00	275,00	304,00	218	

Edit by Technical Office - G.Pedrotti Verified by Quality Assurance - V. Brusa