



**Product Name Surgifix Zerolatex**

<b>Producer</b>
FRA Production S.p.A.
<b>Expiry Date</b>
6 years
<b>Product Description</b>
Elastic netting made of polyamide and elastomeric yarns for medical use. Available in different sizes. Class I medical device - CND code: M03030103.
<b>Intended use</b>
Elastic tubular netting for fixing compresses, gauzes and medications.
<b>Reference Regulations</b>
European Directives 93/42/EC and 2007/47/EC and successive amendments.
<b>Certification</b>
UNI EN ISO 9001:2008. ISO 13485:2012. CE marking.
<b>Warnings</b>
Do not apply directly in contact with the wound. This product is not sterile. Do not put in contact with flames. Do not swallow. Avoid improper use. Keep out of reach of children. This product is to be worn, not wrapped. This product does not contain latex. If the product is too tight, use a larger size.

**Instructions for use**

Wear this product in order to keep sanitary medications in place, choosing the appropriate size.

**Product storage**

Keep away from heat, humidity and light sources.

**Related enclosed documents**

Declaration of conformity to the reference regulations of each batch. Bill of lading.

**Packaging information**

Available in different lengths, as specified in the list of the available product versions.

**Disposal information**

Disposal must be carried out in compliance with related regulations.

In Italy, this product is disposable as "waste from processed textile fibers".

**Guaranteed performances**

Dimensional characteristics of the product, when stretched, are guaranteed.

This product lets skin breathe and allows maximum freedom of movement, though it keeps the dressing in place. Biocompatibility tests have been carried out on the product, in order to prove its use suitability.

**Composition**

86% POLIAMMIDE 14% ELASTOMERO



PRODUCTION S.p.A.

SP000P21  
"Surgifix Zerolatex" elastic netting

Revision 005  
Date 16/03/17

**Product available versions**

Item code	Item description	Sales pack	Number of packs per export carton
00250001	9 SURGIFIX ZEROLATEX TIPO A-I 25 MT	Mt 25,00	20

**Product technical data**

Weight when stretched (g/m) (due to an applied load of of 200g per counter stitch)			Width (mm) when unstretched			Expansion (mm) due to an applied load of 13Kg			Width (mm) after the expansion			Longitudinal elongation (%) (due to a 200g applied load)			Number of counter stitches
std	min	max	std	min	max	std	min	max	std	min	max	std	min	max	
7,00	6,30	7,70	113,0	101,2	124,3	850,0	765,0	935,0	250,0	225,0	275,0	290,00	275,00	304,00	218

Edit by Technical Office - G.Pedrotti

Verified by Quality Assurance - V. Brusa

Via delle Poste 16 - 14010 Dusino San Michele (Italy) - Tel. +390141979911 - Fax +390141979993

\$IFS10

Page 3 of 3

Printed on : 24/05/17