


Adhesive IV dressing REF 4640 (Vnr. 104010)

Version: 06.05.2013, rev.2

Printed: 11.5.2016

<table border="1"> <thead> <tr> <th>REF</th> <th>Size cm</th> <th>Pad size cm</th> </tr> </thead> <tbody> <tr> <td>4640</td> <td>8 x 6</td> <td>1,3x2,4 cm + 1,5x3 cm loose</td> </tr> </tbody> </table>	REF	Size cm	Pad size cm	4640	8 x 6	1,3x2,4 cm + 1,5x3 cm loose	<p>4640</p> 							
REF	Size cm	Pad size cm												
4640	8 x 6	1,3x2,4 cm + 1,5x3 cm loose												
<p>CE0434 The product complies with legal requirements for medical devices</p> <p>EO sterilized according to valid version of ISO 11135 standard</p> <p>NATURAL LATEX FREE</p>	<p>Raw material water repellent, adhesive, spunlace nonwoven 100 % polyester pressure sensitive glue, polyacrylic, latex free viscose/polyester coated paper</p>													
<p>REACH regulations:</p>	<p>The product covered by this data sheet does not contain any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).</p>													
<p>Primary use:</p>	<p>Adhesive dressing for secure fixation of intravascular catheters</p>													
<p>Shelf life:</p>	<p>3 years</p>													
<p>Pack counts & dispenser/carton sizes</p>	<table border="1"> <thead> <tr> <th>REF</th> <th>pcs/ dispenser</th> <th>pcs/ carton</th> <th>box size/ cm WxLxH</th> <th>ctn size/ cm WxLxH</th> </tr> </thead> <tbody> <tr> <td>4640</td> <td>50</td> <td>2400</td> <td>9,4x9,5x15</td> <td>60x40x32</td> </tr> </tbody> </table>				REF	pcs/ dispenser	pcs/ carton	box size/ cm WxLxH	ctn size/ cm WxLxH	4640	50	2400	9,4x9,5x15	60x40x32
REF	pcs/ dispenser	pcs/ carton	box size/ cm WxLxH	ctn size/ cm WxLxH										
4640	50	2400	9,4x9,5x15	60x40x32										
<p>Sales unit:</p>	<p>Transport carton</p>													
<p>Dispenser box</p>	<p>ECF (or TCF) cellulose</p>													
<p>Transport carton:</p>	<p>ECF (or TCF) cellulose</p>													
<p>Packing tape:</p>	<p>PVC-free</p>													
<p>Compliance with legal requirements:</p>	<p>Medical Devices Directive 93/42/EEC, Class Is CE434</p>													
<p>Product safety:</p>	<p>BS EN ISO 10993-parts 5, 10 • cytotoxicity test, PASS • contact sensitization test, PASS • skin irritation test, PASS</p>													
<p>Product performance:</p>	<p>MVTR test –determination of moisture vapor transmission rate - nominal value: > 500 g/m2 during 24 hours</p>													
<p>Quality standards:</p>	<p>The product is manufactured in compliance with ISO13485</p>													
<p>Label information:</p>	<p>in accordance with EN ISO 15223-1, EN 980 and EN 1041</p>													
<p>Storage recommendation:</p>	<p>Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight , other intensive light sources and ozone</p>													
<p>Disposal recommendation</p>	<p>Follow local instructions</p>													

Product Data Sheet

Wound care



Adhesive IV dressing REF 4640 (Vnr. 104010)

Page 2 (2)

Version: 06.05.2013, rev.2

Printed: 11.5.2016

Barcodes: Unit pack	Dispenser	Transport carton
6438129113969 4640 single  6 438129 113969	6438129313963 4640 dispenser  6 438129 313963	6438129513967 4640 carton  6 438129 513967