



Product description

Order number	77.575
Product description	Urine container 3.0 l with transparent inspection strip and graduation, enclosed label with instructions for use

Product characteristics

Label/ Print	With paper label
Graduation	yes
Type of cap	Screw cap
Light protection	yes
Application	Urine collection
Viewing window	yes

Size

Sample volume	3000 ml
Max. volume	3 l
Diameter of opening	83.5 mm
Width of product	130 mm
Length of product	125 mm
Length including cap	246 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Length excluding cap 243 mm

Material & colours

Product material	Polyethylene (PE)
Colour of product	brown
Cap material	Polypropylene (PP)
Colour of cap	yellow

Purity & certification

Product in accordance with MDD or IVDD	IVD
CE certified	CE
Pyrogen-free	no
Safety data sheet	nein
Lot or serial number	Lot number

Packaging

Minimum order qty.	30
Pcs. / inner box	30
Pcs. / box	30
Pcs. / pallet	270
Depth of case	643 mm
Width of case	403 mm
Height of case	515 mm
Case volume	0,1334 cbm
Weight of case	4,21 kg
EAN of inner box	4038917082600
EAN of outer box	4038917029438

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

