



### Product description

Order number	10.252.020
Product description	Urine Monovette®, Luer, 10 ml, 102x15 mm, transparent label, incl. collection tip, 100/inner box individually wrapped sterile

### Product characteristics

Type of collection	Aspiration
Type of preparation	Without preparation
Label/ Print	With plastic label
Colour of print/label	yellow
Graduation	yes
Cap	Luer screw cap
Type of cap	Screw cap
Light protection	no
Application	Urine tube
Base shape	Round base
Volume of work	10 ml
Recommended centrifugation	400 x g
Centrifugation time	5 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Size**

Diameter	15 mm
Length including cap	117 mm
Length excluding cap	102 mm
Length including cap and plunger	127 mm

### **Material & colours**

Colour of product	transparent
Cap material	High Density Polyethylene (HD-PE)
Colour of cap	yellow
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

### **Purity & certification**

Satisfies the requirement	IATA, ADR
Product in accordance with MDD or IVDD	IVD
CE certified	CE
Sterilisation	yes
Pyrogen-free	no
Lot or serial number	Lot number

### **Packaging**

Minimum order qty.	500
Type of inner packaging	individually wrapped and sterile
Pcs. / inner box	100
Pcs. / box	500
Pcs. / pallet	13500
Depth of box	351 mm
Width of box	151 mm
Height of box	171 mm
Depth of case	778 mm
Width of case	368 mm
Height of case	195 mm
Case volume	0,0558 cbm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Weight of case	5,82 kg
EAN of inner box	4038917104432
EAN of outer box	4038917009874

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

