



Product description

Order number	10.252
Product description	Urine Monovette® , Luer, 10 ml, 102x15 mm, transparent label, incl. collection tip, 64/bag

Product characteristics

Type of collection	Aspiration
Type of preparation	Without preparation
Label/ Print	With plastic label
Colour of print/label	yellow
Graduation	yes
Cap	Luer screw cap
Type of cap	Screw cap
Light protection	no
Application	Urine tube
Base shape	Round base
Volume of work	10 ml
Recommended centrifugation	400 x g
Centrifugation time	5 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Size

Diameter	15 mm
Length including cap	117 mm
Length excluding cap	102 mm
Length including cap and plunger	127 mm

Material & colours

Colour of product	transparent
Cap material	High Density Polyethylene (HD-PE)
Colour of cap	yellow
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

Purity & certification

Satisfies the requirement	IATA, ADR
Product in accordance with MDD or IVDD	IVD
CE certified	CE
Pyrogen-free	no
Lot or serial number	Lot number

Packaging

Minimum order qty.	512
Type of inner packaging	Bag
Pcs. / inner box	64
Pcs. / box	512
Pcs. / pallet	24576
Depth of case	388 mm
Width of case	378 mm
Height of case	215 mm
Case volume	0,0315 cbm
Weight of case	4,85 kg
EAN of inner box	4038917279901
EAN of outer box	4038917279772

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

