



Product Name Unislip

Producer
Fra Production S.p.A.

Expiry Date
6 years

Product Description
Elastic net panties made of polyamide, polyester and elastomeric yarns. Class I medical device - CND code: T04010104.

Intended use
Elastic panties intended for medication and incontinence.

Reference Regulations
European Directives 93/42/EC and 2007/47/EC and successive amendments.

Certification
UNI EN ISO 9001:2015. ISO 13485:2016. CE MARKING

Warnings
Do not apply directly in contact with the wound. If the product is too tight, use the next larger size. Do not put in contact with flames. Do not swallow. Avoid improper use. This product is latex-free.

Instructions for use

Wear this product in order to keep sanitary pads or medications in place, choosing the appropriate size. This product is washable many times at 60°C.
Washing instructions: wet clean. Do not chlorine bleach. Do not iron. Do not dry clean; do not remove stains with solvents. Do not tumble dry.

Product storage

Keep away from heat, humidity and light sources.

Related enclosed documents

Declaration of conformity to the reference regulations of each batch. Bill of lading.

Packaging information

The product is sold in bulk or packed inside bags. The specific amounts are written in the "Product available versions" section.

Disposal information

Disposal must be carried out in compliance with related regulations.
In Italy, this product is disposable as "waste from processed textile fibers".

Guaranteed performances

Dimensional characteristics are guaranteed. This product is latex-free.

Composition

PA87% EL12% PES1%



PRODUCTION S.p.A.

SP000R02
"Unislip" elastic panties

Revision 012
Date 5/09/18

Product available versions															
Item code		Item description						Sales pack				Number of bags per carton b			
00560198		UNISLIP S VERDE						PZ. 100/SC				N. 100,00		10	

Product technical data													
Height (cm) when unstretched		Height (cm) when stretched		Length (cm) when unstretched		Length (cm) when stretched due to an applied load of 25Kg		Leg length (cm) when stretched due to an applied load of 25Kg		Crotch length(cm)when stretched due to an applied load of 25Kg		Average liquid containment (g)	Weight (g)
min	max	min	max	min	max	min	max	min	max	min	max		
15,5	17,5	60,0	63,0	30,0	32,0	91,2	96,2	36,6	38,3	18,0	19,5	1060	13,1

Product available versions															
Item code		Item description						Sales pack				Number of bags per carton b			
00560196		UNISLIP M AZZURRO						PZ. 100/SC				N. 100,00		10	

Product technical data													
Height (cm) when unstretched		Height (cm) when stretched		Length (cm) when unstretched		Length (cm) when stretched due to an applied load of 25Kg		Leg length (cm) when stretched due to an applied load of 25Kg		Crotch length(cm)when stretched due to an applied load of 25Kg		Average liquid containment (g)	Weight (g)
min	max	min	max	min	max	min	max	min	max	min	max		
15,5	17,5	60,0	63,0	30,0	32,0	93,2	98,2	37,6	39,3	18,0	19,5	1025	13,4

Product available versions															
Item code		Item description						Sales pack				Number of bags per carton b			
00560197		UNISLIP L MARRONE						PZ. 100/SC				N. 100,00		10	

Product technical data													
Height (cm) when unstretched		Height (cm) when stretched		Length (cm) when unstretched		Length (cm) when stretched due to an applied load of 25Kg		Leg length (cm) when stretched due to an applied load of 25Kg		Crotch length(cm)when stretched due to an applied load of 25Kg		Average liquid containment (g)	Weight (g)
min	max	min	max	min	max	min	max	min	max	min	max		
15,5	17,5	60,0	63,0	32,0	34,0	102,0	107,5	42,0	44,0	18,0	19,5	1400	14,8

Product available versions															
Item code		Item description						Sales pack				Number of bags per carton b			
00560199		UNISLIP XL ARANCIO						PZ. 100/SC				N. 100,00		8	

Product technical data													
Height (cm) when unstretched		Height (cm) when stretched		Length (cm) when unstretched		Length (cm) when stretched due to an applied load of 25Kg		Leg length (cm) when stretched due to an applied load of 25Kg		Crotch length(cm)when stretched due to an applied load of 25Kg		Average liquid containment (g)	Weight (g)
min	max	min	max	min	max	min	max	min	max	min	max		
16,5	18,5	64,0	69,0	33,0	35,0	110,2	115,2	45,1	46,6	20,0	22,0	1200	17,1

Edit by Technical Office - G.Pedrotti Verified by Quality Assurance - V. Brusa

Via delle Poste 16 - 14010 Dusino San Michele (Italy) - Tel. +390141979911 - Fax +390141979993