

# Qufora® Bed Bags | Drainable

Free-flow T-Valve or Lever Valve | Single Packed

Hygienic. Easy-to-use.

The **Qufora® drainable bed bags** can be attached to either a urinary catheter or a urisheath. They are designed to provide additional drainage capacity overnight or for people who are permanently bedridden.

## PRODUCT BENEFITS

- ▶ The Qufora® bed bags have a universal **conical inlet connector** (X-mas tree or step connector) making it easy to attach a urinary catheter or urisheath. The **protection cap** prevents contamination before connecting the tube to the catheter or urisheath. In addition, the cap prevents urine from dropping from the inlet tube during the exchange of the bag.
- ▶ The **flexible connection tubes** provide for safe drainage and allow the bag to be placed below the bladder level.
- ▶ A **non-return valve** ensures that urine will not flow back into the tube, preventing the risk of infection.
- ▶ The Qufora® bed bag (REF 13851202) comes with a fixed and easy-to clean **sampling port** which allows for **needle-free** inspection and eliminates the risk of needle stick injuries. Additionally, there is an **integrated bacteria-proof** and **liquid-tight ventilation** filter for optimal pressure balance. A separate needle-free sampling port connector/protection cap and the Qufora® **hanger** are also included.
- ▶ The **T-Valves** and **Lever valves** have **free-flow** outlets ensuring that clots can freely pass.
- ▶ The Qufora® bed bags can be emptied in a **hygienic** and **safe** way.
- ▶ The front has a **scale** for quick estimation of the urine volume.
- ▶ The Qufora® bed bags have **holes** for the Qufora® hangers or Qufora® floor stands.
- ▶ The Qufora® bed bags are **single packed** with the REF number, lot number and expiration year/month printed on the protective bag.

## Qufora® Accessories (not included):

- ▶ The Qufora® **hangers** are durable and fit most beds.
- ▶ The Qufora® **floor stands** come with four small hooks, making it easy to hang the bed bag.

13853402



Connector with needle-free sample port




T-Valve



Lever valve



Scale



# Qufora® Bed Bags | Drainable

Free-flow T-Valve or Lever Valve | Single Packed

## Product description

- ▶ Qufora® bed bags are available in several versions with different bag volumes or lengths of the inlet tube.
- ▶ The Qufora® bed bags are single-use medical devices.  
**All materials are 100% latex- and DEHP-free.**

## QUFORA® NON-STERILE BED BAGS FEATURES

- Inlet tube: Inner diameter for all bags 4.5mm except 13853402/13853403: 6.6mm and 13851202: 7.6mm
- Non-return valve
- Conical inlet connector (X-mas tree or step)
- T-Valve with white "stopper" (except 13853401)
- Lever valve on 13851202
- Free-flow outlet tube: inner diameter 9.7mm (13853101/13853403: 8.7mm)
- Translucent and white films
- Scale on the front
- Single packed
- Holes for Qufora® hanger or Qufora® floor stand
- 13851202: refer to details on page 3


## MATERIALS

- |   |   |
|---|---|
| • Films                                 | PVC (plasticizer: DOTP/DEHT)                    |
| • Inlet tube                            | PVC (plasticizer: DOTP/DEHT)                    |
| • Needle-free sample port and connector | PVC; Port: TPE                                  |
| • Connector                             | PVC   |
| • Connector cap                         | PE  |
| • Non-return valve                      | PVC (plasticizer: DOTP/DEHT)                    |
| • T-Valve/inlet tube                    | PP/PE; inlet tube: PVC (plasticizer: DOTP/DEHT) |
| • Lever valve                           | LDPE/ABS  |
| • Hanger                                | PP  |
| • Single bag (product)                  | LDPE  |
| • Single bag (hanger)                   | LDPE  |
| • Poly bag (outer)                      | PE  |

## PACKAGING

- |               |   |
|---------------|---|
| • Retail box  | Cardboard   |
| • Shipper box | Corrugated cardboard  |
| • IFU         | One (1) per retail box  |
| • Shelf-life  | Five (5) years  |
| • Storage     | Dry in temperatures between 0°C and 35°C; protected from direct sunlight; humidity level max 80%; not to be exposed to heavy pressure |



**QUFORA® BED BAGS - DRAINABLE  
NON-STERILE - T-VALVE - SINGLE PACKED**

REF	Volume	Tube length	Free-flow outlet tube <sup>3</sup>	White "stopper"	Quantity retail box	Quantity shipper box	GTIN-14 <sup>4</sup>
13753101	1.500 ml	90 cm	Yes	Yes	25	100	15710187004600
13850301	2.000 ml	130 cm	Yes	Yes	10	100	15710187001623
13853101	2.000 ml	90 cm	Yes	Yes	25	100	15710187001562
13853401	2.000 ml	110 cm	Yes	No	25	100	15710187001593
13853402	2.000 ml	110 cm <sup>2</sup>	Yes	Yes	-	100	15710187003634 <sup>5</sup>
13853403 <sup>1</sup>	2.000 ml	110 cm <sup>2</sup>	Yes	Yes	-	100	15710187007311 <sup>5</sup>

All products with X-mas connector (except 13853401: Step connector)

<sup>1</sup>Fixed needle-free sample port connector (X-mas tree connector); <sup>2</sup>Wide inlet tube (X-mas tree connector): Inner diameter: 6.6 mm;

<sup>3</sup>Inner diameter: 9.7mm; <sup>4</sup>GTIN barcode number for retail box; <sup>5</sup>GTIN barcode number for shipper box

**QUFORA® BED BAGS - DRAINABLE  
NON-STERILE - LEVER VALVE - ADJUSTABLE TUBE LENGTH - NEEDLE-FREE SAMPLE PORT - SINGLE PACKED**

REF	Volume	Tube length	Free-flow outlet tube	Separate sample port	Quantity poly bag	Quantity shipper box	GTIN-14 <sup>5</sup>
13851202 <sup>1</sup>	2.000 ml	130 cm <sup>2</sup>	Yes <sup>3</sup>	Yes <sup>4</sup>	10	80	15710187013558

Needle-free sample port, ventilation filter, sheet clip, sliding clamp, hanger and holes for the Qufora® hangers or Qufora® floor stands.

<sup>1</sup>Dog-house for the outlet valve; <sup>2</sup>Wide inlet tube (X-mas tree connector): Inner diameter: 7.6 mm; adjustable tube length; <sup>3</sup>Inner diameter: 8.7mm; <sup>4</sup>Extra/separate needle-free sample port connector (X-mas connector) and protection cap; <sup>5</sup>GTIN barcode number for poly bag

**PRODUCT CLASSIFICATION**

- Medical Devices conforming to MDR 2017/745 in Class I

**APPLIED STANDARDS**

- ISO 13485 Quality Management System
- ISO 8669-2 Urine collection bags  
- Part 2: Requirements and test methods
- GMDN 58922

## Qufora® Accessories (not included)

### QUFORA® BED HANGER



REF	Description	Single packed	Quantity retail box	Quantity shipper box	GTIN-14 <sup>1</sup>
21100101	Bed hanger - fits all beds - easy to operate with one hand	Yes	5	50	15710187001777

### QUFORA® FLOOR STAND



REF	Description	Single packed	Quantity retail box	Quantity shipper box	GTIN-14 <sup>1</sup>
99910003	Floor stand	Yes	1	20	15710187013589

<sup>1</sup>GTIN barcode number for retail box

### PRODUCT CLASSIFICATION

- Medical Devices conforming to MDR 2017/745 in Class I

### APPLIED STANDARDS

- ISO 13485 Quality Management System
- GMDN Bed hanger: 41052; Floor stand: 41052

#### Importer

MBH-International A/S  
Gydevang 28-30  
DK-3450 Allerød  
Denmark

Tel: +45 88 30 05 92  
info@qufora.com  
www.qufora.com