



Product description

Order number 95.1007
Product description Tourniquet

Product characteristics

Colour of print/label white
Application Venous blood collection

Material & colours

Product material Polyamide (PA)
Colour of product red

Purity & certification

Product in accordance with MDD or IVDD Medical product
CE certified CE
Latex-free no
Pyrogen-free no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

Minimum order qty.	1
Type of inner packaging	individually in the bag
Pcs. / inner box	1
Pcs. / box	200
Pcs. / pallet	12800
Depth of case	330 mm
Width of case	340 mm
Height of case	230 mm
Case volume	0,0254 cbm
Weight of case	7 kg
EAN of inner box	4038917084963
EAN of outer box	4038917003728

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

