

BSB X-Ray Cassette Cover

BedsideBag is used under X-ray examinations. The BSB bag with the integrated cushion is pulled over the cassette, the flap is closed and the low-friction bag with the cassette is moved under the patient and positioned in relation to the X-ray tube.

BSB comes with a unique combination of advantages:

- Ease of use and personell efficiency:

BSB enables personell to handle patients quicker, easier and with less effort.

- Patient safety and hygienic standards:

BSB is a hygienic solution which reduces the risk of contamination.

- Patient comfort and relaxation:


BSBs cushioning is relaxing. Relaxed patients are easier to examine.

- Enviromentally friendly:

BSBs material is easy to recycle.

- Reduced stress on the equipment:

BSB is protecting your expensive equipment.



BSB is a win-win solution for both personell and patients.



BSB X-Ray Cassette Cover

BedsideBag

- reduces the work-load for the personell positioning the patient
- protects cassettes and DR-plates from exposure to bodily fluids
- is pressure-reducing, comfortable and insulating
- is made of strong, enviromentally-friendly PE films
- is a proven product used for over 15 years by radiologists i Scandinavia

BedsideBag is an enviromentally-friendly high-quality product.

BSB 2430:

For cassettes with an image format of
Landscape format - Sold in boxes of 50 pcs.

24x30cm / 10x12"

Order no. BSB 2430

BSB 4335:

For cassettes with an image format of
Landscape format - Sold in boxes of 50 pcs.

35x43cm / 14x17"

Order no. BSB 4335

BSB 6353:

For detektorplates with an image format of
Max. dimensions of the detector plate:
Portrait format - Sold in boxes of 50 pcs.

35x43cm / 14x17" and 41x41cm

600(H) x 491(W) x 30(T)mm

Order no. BSB 6353



Copyright CARUM AS 2016

Changes and additions may occur without further notice. For detailed Information on delivery options in your region please contact us.