



STERIKING®

**Sterile Barrier Systems & Products
for Hospitals, Dental Practices & Laundries**


WIPAK

We care that you pack safely!

S - Flat Pouches

For packing individual instruments, smaller sets and other items.*

| Art.code | 2P | 3P | Size (mm) | Pieces/case |
|----------|----|----|-----------|-------------|
| S 17 | • | | 50 x 200 | 3.600 |
| S 1 | • | • | 50 x 250 | 2.400 |
| S 24 | • | • | 75 x 150 | 2.400 |
| S 2 | • | • | 75 x 200 | 2.400 |
| S 23 | • | | 75 x 230 | 1.800 |
| S 4 | • | • | 75 x 270 | 3.600 |
| S 18 | • | • | 75 x 300 | 4.200 |
| S 22 | • | | 75 x 520 | 2.400 |
| S 25 | • | | 100 x 150 | 2.400 |
| S 3 | • | • | 100 x 200 | 1.800 |
| S 5 | • | • | 100 x 270 | 3.000 |
| S 8 | • | • | 100 x 300 | 1.200 |
| S 19 | • | | 100 x 350 | 1.200 |
| S 9 | • | • | 100 x 400 | 1.200 |
| S 12 | • | | 100 x 570 | 1.800 |
| S 27 | • | | 120 x 400 | 1.200 |
| S 15 | • | • | 150 x 200 | 1.200 |
| S 6 | • | • | 150 x 270 | 1.800 |
| S 20 | • | • | 150 x 300 | 1.800 |
| S 26 | • | • | 150 x 350 | 1.800 |
| S 13 | • | • | 150 x 400 | 1.200 |
| S 28 | • | | 150 x 520 | 1.200 |
| S 34 | • | | 160 x 600 | 1.200 |
| S 33 | • | | 180 x 320 | 1.000 |
| S 7 | • | • | 205 x 270 | 600 |
| S 7 A | • | | 205 x 340 | 1000 |
| S 10 | • | • | 205 x 400 | 600 |
| S 10 A | • | | 205 x 440 | 600 |
| S 35 | • | | 205 x 600 | 1.200 |
| S 2070 | • | | 205 x 700 | 600 |
| S 16 | • | • | 250 x 380 | 1.200 |
| S 11 | • | • | 250 x 500 | 600 |
| S 36 | • | | 250 x 600 | 500 |
| S 37 | • | | 250 x 700 | 600 |
| S 29 | • | | 270 x 350 | 600 |
| S 30 | • | • | 270 x 440 | 600 |
| S 21 | • | • | 300 x 500 | 600 |
| S 14 | • | • | 300 x 570 | 600 |
| S 38 | • | | 320 x 500 | 600 |
| S 31 | • | | 420 x 500 | 500 |
| S 32 | • | • | 420 x 600 | 500 |



B - Pouches with Gusset

For more voluminous items.*

| Art.code | 2P | Size (mm) | Pieces/case |
|----------|----|-----------------|-------------|
| B 35 | • | 100 x 50 x 300 | 500 |
| B 30 | • | 100 x 50 x 360 | 500 |
| B 31 | • | 150 x 50 x 400 | 500 |
| B 32 | • | 150 x 50 x 460 | 500 |
| B 36 | • | 200 x 55 x 400 | 250 |
| B 33 | • | 200 x 55 x 500 | 250 |
| B 37 | • | 250 x 65 x 480 | 250 |
| B 34 | • | 300 x 80 x 550 | 250 |
| B 38 | • | 420 x 100 x 570 | 250 |

R - Flat Rolls *

| Art.code | 2P | 3P | Size (mm) | Pieces/case |
|----------|----|----|-------------|-------------|
| R 39 | • | • | 50 : 200 m | 5 |
| R 40 | • | • | 75 : 200 m | 4 |
| R 41 | • | • | 100 : 200 m | 3 |
| R 41A | • | • | 120 : 200 m | 2 |
| R 42 | • | • | 150 : 200 m | 2 |
| R 43 | • | • | 200 : 200 m | 1 |
| R 44 | • | • | 250 : 200 m | 1 |
| R 45 | • | • | 300 : 200 m | 1 |
| R 46 | • | • | 350 : 200 m | 1 |
| R 47 | • | • | 400 : 200 m | 1 |

| Art.code | 2P | 3P | Size (mm) | Pieces/case |
|----------|----|----|-------------|-------------|
| R 39/100 | • | | 50 : 100 m | 8 |
| R 40/100 | • | | 75 : 100 m | 8 |
| R 41/100 | • | | 100 : 100 m | 4 |
| R 42/100 | • | | 150 : 100 m | 2 |
| R 43/100 | • | | 200 : 100 m | 2 |
| R 44/100 | • | | 250 : 100 m | 2 |
| R 45/100 | • | | 300 : 100 m | 1 |
| R 47/100 | • | | 400 : 100 m | 1 |
| R 49/100 | • | | 500 : 100 m | 2 |

| Art.code | 2P | 3P | Size (mm) | Pieces/case |
|----------|----|----|------------|-------------|
| R 39/31 | • | | 50 : 31 m | 10 |
| R 40/31 | • | | 75 : 31 m | 10 |
| R 41/31 | • | | 100 : 31 m | 10 |
| R 42/31 | • | | 150 : 31 m | 10 |
| R 43/31 | • | | 200 : 31 m | 6 |
| R 44/31 | • | | 250 : 31 m | 6 |
| R 45/31 | • | | 300 : 31 m | 6 |

RB - Rolls with Gusset *

| Art.code | 2P | 3P | Size (mm) | Pieces/case |
|----------|----|----|------------------|-------------|
| RB 50 | • | • | 75 x 25 : 100 m | 4 |
| RB 51 | • | • | 100 x 50 : 100 m | 3 |
| RB 52 | • | • | 150 x 50 : 100 m | 2 |
| RB 53 | • | • | 200 x 55 : 100 m | 1 |
| RB 54 | • | • | 250 x 65 : 100 m | 1 |
| RB 55 | • | • | 300 x 80 : 100 m | 1 |
| RB 56 | • | • | 350 x 80 : 100 m | 1 |
| RB 57 | • | • | 400 x 80 : 100 m | 1 |

* 2P = steam and EO gas, 3P = steam, EO and FO.

SMX - High Performance ProWraps

SMX-1 Blue Light

For small and medium-size trays and sets.

| Art.code | Size (mm) | Pieces/case |
|----------|-------------|-------------|
| SMX 1060 | 600 x 600 | 480 |
| SMX 1075 | 750 x 750 | 200 |
| SMX 1090 | 900 x 900 | 160 |
| SMX 1100 | 1000 x 1000 | 160 |
| SMX 1120 | 1200 x 1200 | 120 |

SMX-2 Blue Regular

For medium and large-size trays and sets and for container inner wraps.

| Art.code | Size (mm) | Pieces/case |
|----------|-------------|-------------|
| SMX 2060 | 600 x 600 | 480 |
| SMX 2075 | 750 x 750 | 200 |
| SMX 2090 | 900 x 900 | 120 |
| SMX 2100 | 1000 x 1000 | 120 |
| SMX 2120 | 1200 x 1200 | 120 |

SMX-2+ Green Regular

Optimal choice for inner wrap, as it offers a reliable sterile field.

| Art.code | Size (mm) | Pieces/case |
|-----------|-------------|-------------|
| SMX 2+060 | 600 x 600 | 480 |
| SMX 2+075 | 750 x 750 | 200 |
| SMX 2+090 | 900 x 900 | 120 |
| SMX 2+100 | 1000 x 1000 | 120 |
| SMX 2+120 | 1200 x 1200 | 120 |

SMX-3 Blue Heavy Duty

Heavier structure for the most demanding applications.

| Art.code | Size (mm) | Pieces/case |
|----------|-------------|-------------|
| SMX 3075 | 750 x 750 | 120 |
| SMX 3090 | 900 x 900 | 120 |
| SMX 3100 | 1000 x 1000 | 120 |
| SMX 3120 | 1200 x 1200 | 104 |

SMX-4 Blue Super Duty

Especially high thickness for orthopedic sets and other heavy-duty uses.

| Art.code | Size (mm) | Pieces/case |
|----------|-------------|-------------|
| SMX 4075 | 750 x 750 | 120 |
| SMX 4090 | 900 x 900 | 104 |
| SMX 4100 | 1000 x 1000 | 104 |
| SMX 4120 | 1200 x 1200 | 80 |

NWB - Nonwoven Soft Blue

For general use, for laundry and instrument sets.

| Art.code | Size (mm) | Pieces/case |
|----------|-------------|-------------|
| NWB 50 | 500 x 500 | 500 |
| NWB 60 | 600 x 600 | 500 |
| NWB 75 | 750 x 750 | 250 |
| NWB 90 | 900 x 900 | 250 |
| NWB 100 | 1000 x 1000 | 250 |
| NWB 120 | 1200 x 1200 | 100 |
| NWB 137 | 1370 x 1370 | 100 |

NWG - Nonwoven Super Soft Green

Superior softness and toughness for demanding applications.

| Art.code | Size (mm) | Pieces/case |
|--|----------------------|-------------|
| NWG 50 | 500 x 500 | 400 |
| NWG 60 | 600 x 600 | 400 |
| NWG 75 | 750 x 750 | 200 |
| NWG 90 | 900 x 900 | 150 |
| NWG 100 | 1000 x 1000 | 150 |
| NWG 120 | 1200 x 1200 | 100 |
| Available as interleaved with SMX Wrap-grade | | |
| NWGL | Nonwoven Green Light | |



SPC - Sterilization Paper Crepe

Special grade paper creped to give increased flexibility.

| Art.code | Size (mm) | Pieces/case |
|-------------|-------------|-------------|
| SPC 40 | 400 x 400 | 500 |
| SPC 50 | 500 x 500 | 500 |
| SPC 60 | 600 x 600 | 500 |
| SPC 75 | 750 x 750 | 250 |
| SPC 90 | 900 x 900 | 250 |
| SPC 100 | 1000 x 1000 | 250 |
| SPC 120 | 1200 x 1200 | 100 |
| SPC 100/140 | 1000 x 1400 | 100 |

SCB - Super Crepe Blue

Soft crepe - for improved drapeability.

| Art.code | Size (mm) | Pieces/case |
|----------|-------------|-------------|
| SCB 50 | 500 x 500 | 500 |
| SCB 60 | 600 x 600 | 500 |
| SCB 75 | 750 x 750 | 250 |
| SCB 90 | 900 x 900 | 250 |
| SCB 100 | 1000 x 1000 | 250 |
| SCB 120 | 1200 x 1200 | 100 |

Absorbent Trayliners

For use in surgical instrument trays to absorb excessive condensation.

| TL-Standard | Size (mm) | Pieces/case |
|---------------------------|-----------|-------------|
| TL 25 | 250 x 300 | 1000 |
| TL 40 | 400 x 500 | 500 |
| TL 50 | 500 x 300 | 500 |
| TLWE-Trayliners with edge | | |
| TLWE 2525 | 250 x 250 | 500 |
| TLWE 4825 | 480 x 250 | 500 |
| TLWE 5030 | 500 x 300 | 500 |
| TLWE-Trayliners nonwoven | | |
| TLNW 1525 | 150 x 250 | 3.200 |
| TLNW 3025 | 300 x 250 | 1.200 |
| TLNW 3040 | 300 x 400 | 1.200 |
| TLNW 3050 | 300 x 500 | 600 |



Indicator Strips

AC = steam, EO-gas, FO-gas and
PL = low temperature plasma.

NEW!

| Art.code | Strips/dispenser | Pieces/case |
|------------|------------------|-------------|
| AC 6 | 250 | |
| AC 250/500 | 250 | 1.000 |
| EO 250 | 250 | 1.000 |
| FO 250/500 | 250 | |
| PL 500 | | 500 |

Daily Test Sheet

For Bowie & Dick Test with a cloth package.

| Art.code | Pieces/case |
|----------|-------------|
| DT 50 | 50 |

Daily Test Pack

The ready made test packs for Bowie & Dick Test.

| Art.code | Pieces/case |
|----------|-------------|
| DCC2 | 20 |
| DC 20 | 20 |

Helix Challenge Test

For a charge control and as a steam penetration test for hollow instruments with small lumina.

| Art.code | Strips/dispenser | Pieces/case |
|----------|------------------|-------------|
| HX 100 | 100 | 25 |
| HX 400 | 400 | 25 |

Tapes

Indicator Tapes, S=steam, E=EO-gas and H=hot-air indicator

| Art.code | Width/length | Pieces/case |
|-----------------------|--------------|-------------|
| ITS 12 G | 12 mm x 50 m | 72 |
| ITS 19 G | 19 mm x 50 m | 48 |
| ITS 25 G | 25 mm x 50 m | 36 |
| JTB 19 S G (B = blue) | 19 mm x 50 m | 48 |
| JTB 25 S G (B = blue) | 25 mm x 50 m | 36 |
| ITE 19 | 19 mm x 50 m | 48 |
| ITH 19 | 19 mm x 50 m | 48 |

Closing Tapes with no imprint

| Art.code | Width/length | Pieces/case |
|---------------------|--------------|-------------|
| ITN 19 | 19 mm x 50 m | 48 |
| ITN 25 | 25 mm x 50 m | 36 |
| JTB 19 N (B = blue) | 19 mm x 50 m | 48 |
| JTB 25 N (B = blue) | 25 mm x 50 m | 36 |

Plasma Indicator Tape Strip of DuPont™ Tyvek®

| Art.code | Width/length | Pieces/case |
|----------|----------------|-------------|
| ITP 20 | 20 mm x 150 mm | 200 |

Tape dispenser for one roll

| Art.code | Pieces/case |
|----------|-------------|
| TDS 25 | 1 |

SS - Standard Self Seal Pouches

with adhesive closure. 2P = steam and EO gas. **NEW UNITS!**

| Art.code | Size (mm) | Pieces/case | Pieces/case |
|----------|-----------|-------------|-------------|
| SS 1 | 90 x 200 | 1.200 | 1.000 |
| SS 2 | 90 x 250 | 1.200 | 1.000 |
| SS 3 | 90 x 570 | 1.200 | 1.000 |
| SS 4 | 130 x 270 | 1.200 | 1.000 |
| SS 4 A | 130 x 380 | 1.200 | 1.000 |
| SS 5 | 200 x 350 | 1.200 | 1.200 |
| SS 5 A | 190 x 330 | 1.200 | 1000 |
| SS 6 | 250 x 400 | 1.200 | 600 |
| SS 7 | 300 x 450 | 600 | 600 |

SST - Self Seal Plus Pouches

with an attached internal indicator strip.

| Art.code | Size (mm) | Pieces/case |
|----------|-----------|-------------|
| SST 1 | 90 x 200 | 1.000 |
| SST 2 | 90 x 200 | 1.000 |
| SST 3A | 100 x 640 | 1.000 |
| SST 4 | 130 x 270 | 1.000 |
| SST 4A | 130 x 380 | 1.000 |
| SST 5A | 190 x 330 | 1.000 |
| SST 6 | 250 x 400 | 600 |
| SST 7 | 300 x 450 | 600 |

SNW - TEX-Line Pouches

Peel pouches with nonwoven and film for packing heavy and bulky items. 3P = steam, EO and FO

| Art.code | Size (mm) | Pieces/case |
|----------|-----------|-------------|
| SNW 2040 | 205 x 400 | 500 |
| SNW 2744 | 270 x 440 | 500 |
| SNW 3250 | 320 x 500 | 500 |
| SNW 4260 | 420 x 600 | 300 |
| SNW 5065 | 500 x 650 | 200 |
| SNW 5075 | 500 x 750 | 200 |
| SNW 5775 | 570 x 750 | 200 |

CB - Dust Cover Bags

for use after sterilization.

| Art.code | Size (mm) | Pieces/case |
|----------|-----------|-------------|
| CB 1 | 150 x 270 | 1.000 |
| CB 2 | 205 x 300 | 1.000 |
| CB 3 | 250 x 400 | 1.000 |
| CB 4 | 300 x 500 | 500 |
| CB 5 | 400 x 550 | 500 |
| CB 6 | 420 x 700 | 500 |

SSDC - Self Seal Dust Cover Bags

with Adhesive Closure

| Art.code | Size (mm) | Pieces/case |
|----------|-----------|-------------|
| SSDC 10 | 150 x 270 | 500 |
| SSDC 20 | 200 x 350 | 500 |
| SSDC 30 | 250 x 400 | 500 |
| SSDC 40 | 300 x 550 | 200 |
| SSDC 60 | 420 x 700 | 200 |

LTS - Flat Tyvek® Pouches

Made with DuPont™ Tyvek® for low temperature sterilization.

Printed with hydrogen plasma process indicator. **NEW UNITS!**


| Art.code | Size (mm) | Pieces/case | Pieces/case |
|----------|-----------|-------------|-------------|
| LTS 7520 | 75 x 200 | 1.200 | 1.000 |
| LTS 1025 | 100 x 250 | 1.200 | 1.000 |
| LTS 1530 | 150 x 300 | 1.200 | 1.000 |
| LTS 1644 | 160 x 440 | 1.200 | 600 |
| LTS 1660 | 160 x 600 | 600 | 600 |
| LTS 2038 | 205 x 380 | 1.200 | 1.200 |
| LTS 2538 | 250 x 380 | 1.200 | 600 |
| LTS 2550 | 250 x 500 | 600 | 600 |

LTSS - Self Seal Tyvek® Pouches

| Art.code | Size (mm) | Pieces/case |
|----------|-----------|-------------|
| LTSS 1 | 90 x 200 | 1.000 |
| LTSS 2 | 90 x 250 | 1.000 |
| LTSS 4 | 130 x 270 | 1.000 |
| LTSS 4 A | 130 x 380 | 1.000 |
| LTSS 5 A | 190 x 330 | 1.000 |
| LTSS 6 | 250 x 400 | 600 |
| LTSS 7 | 300 x 450 | 600 |

LTR - Flat Tyvek® Rolls

| Art.code | Size (mm) | Pieces/case |
|----------|-----------|-------------|
| LTR 40 | 75 x 100 | 1 |
| LTR 41 | 100 x 100 | 1 |
| LTR 42 | 150 x 100 | 1 |
| LTR 43 | 200 x 100 | 1 |
| LTR 43 A | 225 x 100 | 1 |
| LTR 44 | 250 x 100 | 1 |
| LTR 45 | 300 x 100 | 1 |
| LTR 46 | 350 x 100 | 1 |
| LTR 47 | 400 x 100 | 1 |
| LTR 49 | 500 x 100 | 1 |





Sealing Machines

NEW! Model: RS 120A Rotosealer

Model: RS 220A Advanced Rotosealer with printer

Model: IS 250 Impulse sealer with cutting device

Cutting Devices in 3 different widths

• CD 600, CD 900, CD 1200

Seal Control Sheet

For operational qualification of the sealing process.

| Art.code | Strips/dispenser | Pieces/case |
|-----------------|------------------|-------------|
| SC 250 (paper) | 250 | 1.000 |
| TSC 200 (Tyvek) | 200 | 800 |

NEW!

Non-toxic Markers

Waterproof, heat-resistant markers for sterilization packages. 10/dispenser.

| Art.code | Color Pieces | Size (mm) | Pieces/case |
|----------|--------------|-----------|-------------|
| MAR 1 | ■ | 0.75 | 100 |
| MAR 1T | ■ | 1.00 | 100 |
| MAR 2 | ■ | 0.75 | 100 |
| MAR 3 | ■ | 0.75 | 100 |
| MAR 4 | ■ | 0.75 | 100 |

PB - Paper Bags with Gusset

Heat-sealable lacquer top and steam indicator.

| Art.code | Size (mm) | Pieces/case |
|----------|-----------------|-------------|
| PB 1 | 90 x 50 x 125 | 1.000 |
| PB 2 | 110 x 30 x 190 | 1.000 |
| PB 3 | 125 x 50 x 250 | 1.000 |
| PB 4 | 140 x 75 x 250 | 500 |
| PB 5 | 140 x 50 x 330 | 500 |
| PB 6 | 165 x 50 x 380 | 500 |
| PB 7 | 180 x 95 x 380 | 500 |
| PB 8 | 190 x 65 x 330 | 500 |
| PB 9 | 250 x 100 x 380 | 500 |
| PB 10 | 300 x 75 x 530 | 250 |

HR - Hot Air Rolls

Special polyamide tubing for sterilization by dry heat.

| Art.code | Width/length | Pieces/case |
|----------|----------------|-------------|
| HR 39 | 50 mm x 200 m | 5 |
| HR 40 | 75 mm x 200 m | 4 |
| HR 41 | 100 mm x 200 m | 3 |
| HR 42 | 150 mm x 200 m | 2 |
| HR 43 | 200 mm x 200 m | 1 |
| HR 44 | 250 mm x 200 m | 1 |

Registration Labels

For registering the sterilized items, enabling recall of all items if necessary.

PP-based labels, resistant to disinfection, washing and sterilization

| Art.code | | Pieces/case |
|----------|----------------------|-------------|
| TRLAB | Tray Label unprinted | 1.000 |
| TRLABPR | Tray Label printed | 1.000 |

Paper label with adhesive backing, suitable as a bar code carrier

| Art.code | | Pieces/case |
|----------|-----------|-------------|
| TRTAG | Tray Tag | 1.000 |
| ENV 1 | Envelopes | 1.000 |

STERIKING®

Wipak Oy
P.O. Box 45
Wipaktie 2
FI-15561 Nastola, Finland
Tel: +358-20510 311
Fax: +358-20510 3333
www.wipak.com
www.steriking.info
e-Mail: steriking@wipak.com



WIPAK

We reserve the right to make changes in the products.