



Irrigation Syringes

3-piece syringes for irrigation purposes



Omnifix® Irrigation Syringe

3-piece syringe with catheter tip



1



2

Product Features

- Highly transparent barrel
- High contrast graduation
- Catheter tip (1)
- Integrated plunger backstop (2)
- Graduation extended beyond nominal volume
- Barrel and plunger made of polypropylene
- Plunger stopper with double sealing ring made of polyisoprene
- Not manufactured with PVC, Latex, DEHP and BPA
- Sterile single packaging
- Box of 100 pcs.

User Benefits

- Highly transparent barrel with high contrast graduation for easy reading and exact dosing
- Plunger stopper with double sealing ring allows slow aspiration and injection
- Safe plunger backstop helps to prevent unintentional plunger withdrawal and the loss of any liquids

Size	Cone	Graduation	Total Graduated Capacity	Units Per Box	Code No. (REF)
50 ml	Centric	1 ml	60 ml	100 pcs.	4613503F