

# PRODUCT DATA SHEET

Omnifix® (50 ml), Irrigation Syringe(4613503F)  
3-piece single-use irrigation syringe



Attribute Name	Attribute Value
Brand name	Omnifix®
Product name (HC)	Omnifix®
Product type	Syringe
Application	Irrigation
Construction of syringe	3-piece
Style of syringe cone	Catheter fitting
Position of syringe cone	centric
Nominal capacity	50.0 ml
Total graduated capacity	60.0 ml
Max. usable capacity	60.00 ml
Scale value	1.0000 ml
Light-protected	No
Color of light protection	-
Type of packaging	Blister
Packaging contains PVC	No

Attribute Name	Attribute Value
Packaging contains Latex	No
Packaging contains DEHP	No
Product (unpack.) contains PVC	No
Product (unpack.) cont. Latex	No
Product (unpack.) cont. DEHP	No
Product contains PVC	No
Product contains Latex	No
Product contains DEHP	No
Product contains Nickel	No
Product contains BPA	No
ISO Code material	Polypropylene
Lubricants	Silicon fluid
Reference norm	ISO 11135-1
Sterility	sterile
Type of Sterilization	EO
Type of Sterilizat. (allowed)	EO
Disposable / Reusable	disposable
eCl@ss Description (Rel. 10)	Irrig./bladder syringe (disp.)
eCl@ss Code (Rel. 10)	34220107
eCl@ss Description (Rel. 11)	Irrig./bladder syringe (disp.)
eCl@ss Code (Rel. 11)	34220107
UMDNS-Version	1.0
UMDNS	13942-Syringes, Irrigating