

# Venofix<sup>®</sup> A

Winged IV needle for optimum venipuncture



Vascular Access

# Venofix® A

## Easy and accurate venipuncture



The Venofix® A product range is comprising gauge sizes from 27G to 19G to be suitable for a wide field of indications – short-term infusion, injection, transfusion and blood sampling.



The flexible interlocking wings and the micro siliconized precision bevel allow an appropriate puncture angle for various venous conditions. Soft wings facilitate a safe fixation on the patient's skin.



### Enhanced short-term infusion therapy

The DEHP-free tubing prevents from kinking to avoid an interruption of infusion therapy. Venofix® A is appropriate for modern infusion therapy in connection with pressure infusion equipment, as it meets the requirements of ISO 8536-9.

### Product specifications

|   | Size (mm) / (G) | Useful length (mm) | Flow rate (ml/min) | Priming volume (ml) | Code No. (REF) |
|---|-----------------|--------------------|--------------------|---------------------|----------------|
| ■ | 0.40 / 27       | 10                 | 2                  | 0.38                | 4056388        |
| ■ | 0.50 / 25       | 15                 | 4                  | 0.39                | 4056370        |
| ■ | 0.55 / 24       | 20                 | 8                  | 0.39                | 4056361        |
| ■ | 0.65 / 23       | 20                 | 11                 | 0.39                | 4056353        |
| ■ | 0.70 / 22       | 20                 | 23                 | 0.39                | 4056345        |
| ■ | 0.80 / 21       | 20                 | 25                 | 0.39                | 4056337        |
| ■ | 1.10 / 19       | 20                 | 43                 | 0.41                | 4056310        |

Tubing length: 30 cm, Luer-Lock fitting with closing cone

Sales unit: 200 pieces