

Arm cover
REF 210750(Vnr. 702293)

Page 1 (2)

Version: 09.08.2013, rev 1

Printed 26.5.2014


REF	Size	
210750	21 x 43 cm	 <p style="text-align: center;">210750</p>
Color & features 210750 white, with elasticated cuff <div style="border: 1px solid blue; padding: 5px; display: inline-block; color: blue; font-weight: bold;">[NATURAL LATEX FREE]</div>		
Raw material: 210750		PE film, elastic both round wrist and upper arm
REACH regulations:		The products covered by this data sheet do not contain any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Features:		Disposable arm cover is a practical and economical product for infection control. Also suited for non-medical uses.
Pack count:		210745 2 pcs/bag, 50 pcs/box, 400 pcs/transport carton 210750 10 pcs/roll, 10 rolls/bag, 2000 pcs/transport carton
Sales unit:		Transport carton
Dispenser box		PE bag 300 x 200 mm
Transport carton:		Cardboard 380 x 250 x 310 mm
Compliance with legal requirements:		Personal Protective equipment 89/686/EEC, cat I
Quality standards:		The product is manufactured in compliance with ISO13485
Label information:		In accordance with EN ISO 15223-1, EN 980 and EN 1041
Storage recommendation:		Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone
Disposal recommendation		Arm & packaging material to energy recovery or mixed waste

Arm cover
REF 210750(Vnr. 702293)

Page 2 (2)

Version: 09.08.2013, rev 1

Printed 26.5.2014

Barcode	Dispenser	Transport carton
210750	N/A	6438129552607 210750 carton  6 438129 552607