

3M™ Littmann® Stethoscopes



Repairs made easy!

3M™ Littmann® Stethoscope Spare Parts Kits

Now, stethoscope customization and simple repairs are easier than ever – with our new line of 3M™ Littmann® Stethoscope Spare Parts Kits. These kits are conveniently packaged to correspond to our stethoscope models and meet the same demanding standards – letting clinicians get back to work faster.

Cat #	Product Name	Kit Includes	UPC 007073-8777-
40003	3M™ Littmann® Stethoscope Spare Parts Kit, Cardiology III™, Black, 40003	<ul style="list-style-type: none"> • Snap Tight Soft-Sealing Eartips – Small • Tunable Rim/Diaphragm Assembly – Adult • Rim, Snap-on – Pediatric • Tunable Diaphragm – Pediatric 	2949
40004	3M™ Littmann® Stethoscope Spare Parts Kit, Cardiology III™, Gray, 40004	<ul style="list-style-type: none"> • Snap Tight Soft-Sealing Eartips – Small • Tunable Rim/Diaphragm Assembly – Adult • Rim, Snap-on – Pediatric • Tunable Diaphragm – Pediatric 	2956
40005	3M™ Littmann® Stethoscope Spare Parts Kit, Classic II S.E., Black, 40005	<ul style="list-style-type: none"> • Nonchill Bell Sleeve • Snap Tight Soft-Sealing Eartips – Small • Tunable Rim/Diaphragm Assembly – Adult 	2963
40006	3M™ Littmann® Stethoscope Spare Parts Kit, Classic II S.E., Gray, 40006	<ul style="list-style-type: none"> • Nonchill Bell Sleeve • Snap Tight Soft-Sealing Eartips – Small • Tunable Rim/Diaphragm Assembly – Adult 	2970
40011	3M™ Littmann® Stethoscope Spare Parts Kit, Master Cardiology™, Black, 40011	<ul style="list-style-type: none"> • Snap Tight Soft-Sealing Eartips – Small • Tunable Rim/Diaphragm Assembly 	2987
40012	3M™ Littmann® Stethoscope Spare Parts Kit, Classic II Pediatric Diaphragm Assembly, 40012	<ul style="list-style-type: none"> • Black Snap-on Rim • Floating Diaphragm • Gray Snap-on Rim 	2994
40013	3M™ Littmann® Stethoscope Spare Parts Kit, Classic II Infant Assembly, 40013	<ul style="list-style-type: none"> • Black Snap-on Rim • Floating Diaphragm • Gray Snap-on Rim 	3007
40016	3M™ Littmann® Stethoscope Spare Parts Kit, Classic III™, Black, 40016	<ul style="list-style-type: none"> • Snap Tight Soft-Sealing Eartips – Small • Tunable Single Piece Diaphragm – Adult • Tunable Single Piece Diaphragm – Pediatric 	3014
40017	3M™ Littmann® Stethoscope Spare Parts Kit, Classic III™, Gray, 40017	<ul style="list-style-type: none"> • Snap Tight Soft-Sealing Eartips – Small • Tunable Single Piece Diaphragm – Adult • Tunable Single Piece Diaphragm – Pediatric 	3021
40018	3M™ Littmann® Stethoscope Spare Parts Kit, Master Cardiology™, Gray, 40018	<ul style="list-style-type: none"> • Snap Tight Soft-Sealing Eartips – Small • Tunable Rim/Diaphragm Assembly 	3038
40019	3M™ Littmann® Stethoscope Spare Parts Kit, Cardiology S.T.C., Black, 40019	<ul style="list-style-type: none"> • Snap Tight Soft-Sealing Eartips – Small • Tunable Rim/Diaphragm Assembly 	3045
40020	3M™ Littmann® Stethoscope Spare Parts Kit, Lightweight II S.E., Black, 40020	<ul style="list-style-type: none"> • Snap Tight Soft-Sealing Eartips – Small • Tunable Rim/Diaphragm Assembly • Nonchill Bell Sleeve 	3052
40021	3M™ Littmann® Stethoscope Spare Parts Kit, Lightweight II S.E., Light Brown, 40021	<ul style="list-style-type: none"> • Snap Tight Soft-Sealing Eartips – Small • Tunable Rim/Diaphragm Assembly • Nonchill Bell Sleeve 	3069
40022	3M™ Littmann® Stethoscope Spare Parts Kit, Master Classic™, Black, 40022	<ul style="list-style-type: none"> • Snap Tight Soft-Sealing Eartips – Small • Tunable Rim/Diaphragm Assembly 	3076
40023	3M™ Littmann® Stethoscope Spare Parts Kit, Master Classic™, Gray, 40023	<ul style="list-style-type: none"> • Snap Tight Soft-Sealing Eartips – Small • Tunable Rim/Diaphragm Assembly 	3083



3M Health Care
 3M Infection Prevention Division
 2510 Conway Ave., St. Paul, MN 55144 USA
 1-800-228-3957 | www.littmann.com

Please recycle. Printed in U.S.A.
 © 3M 2015. All rights reserved.
 70-2011-5754-5



3M, Littmann, the L Littmann logo, L, Classic III, Master Cardiology and Cardiology III are trademarks of 3M. Used under license in Canada.

