



8033651EU

Armatec™ microscope drape for Leica™ OH3, OH4, OHS, OH1, MS-2, MS-3, M520/M525, featuring angled lens with click closure

SIZE: **117 x 305 cm**

QUANTITY: **10**

CONTAINS LATEX: **No**

STERILE: **Yes**

STERILIZATION METHOD: **EO**

MATERIAL: **Polyethylene**

GTIN PRIMARY: **50748426081535**

GTIN SECONDARY: **00748426081530**

GTIN SALES: **50748426081535**

COMMODITY CODE: **3926909790**

INTENDED USE: **Equipment Drapes are used as a protective covering for various surgical and/or nonsurgical instruments/equipment [e.g., a surgical or endoscopic tray cover or liner, an operating room (OR) equipment**

SHELF LIFE: **5 Years**

DEHP/PHTALATES: **No**

GMDN CODE: **12535**

CE NOTIFYING BODY: **CE0044**

MEDICAL CLASSIFICATION: **Is**

BOX LENGTH (MM): **450**

BOX WIDTH (MM): **300**

BOX HEIGHT: **285**

BOX M3: **0.038475**

BOX WEIGHT (GRAM): **3600**

COUNTRY OF ORIGIN: **Dominican Republic**

LABEL LANGUAGES: **BG,CS,DA,DE,EL,EN,ES,ET,FI,FR,HR,HU,IS,IT,LT,LV,NL,NO,PL,PT,RO,RU,SK,SL,SR,SV,TR**

2X STERILISATION ALLOWED?: **Yes**

ADDRESS OF LEGAL MANUFACTURER: **1 Ecolab Place, St. Paul, MN 55102 U.S.A.**

BISPHENOL A: **No**

IFU: **Yes**

PVC: **No**

REACH COMPLIANCE: **Yes**

SINGLE USE: **Yes**

CASES PER PALLET CONVERSION: **39**

INFORMATION LAST UPDATED ON: **21 Jun 2022**
