



# USON PLAST A/S

Care for people

Safety containers and medical supplies for the health care sector



## Præparatbeholder - 3500 ml

Gul beholder med gult låg produceret af polypropylen



### VARENØR:

# 902

### DIMENSIONER

Diameter (top): 17 cm  
Højde: 17 cm

### EMBALLERING

Antal pr. karton: 45 stk.  
Antal pr. palle: 900 stk.



ISO 9001  
BUREAU VERITAS  
Certification

USON PLAST A/S  
is certified by  
BVQI according to  
DS/EN ISO 9001.