

PRODUCT DATA SHEET

STERICAN 0,60X60 G23X2 3/8 NEURAL THERAP (4665600)

Single-use hypodermic needles



Attribute Name	Attribute Value
Brand name	Sterican®
Product name (HC)	Sterican®
Product type	Needle
Application	Neutral therapy
Connection	Luer-Lock/female
Type of bevel	long
Outer-ø cannula	0.60 mm
Gauge (cannula)	23
Color	blue
Length of cannula (metric)	60.00 mm
Length of cannula (inch)	2 3/8"
Residual volume	0.061 ml
Type of packaging	Blister

Attribute Name	Attribute Value
Product (unpack.) contains PVC	No
Product (unpack.) cont. DEHP	No
Product (unpack.) cont. Latex	No
Packaging contains PVC	No
Packaging contains DEHP	No
Packaging contains Latex	No
Product contains PVC	No
Product contains DEHP	No
Product contains Latex	No
Product contains Nickel	Yes
Sterility	sterile
Type of Sterilization	EO
Disposable / Reusable	disposable