

PRODUCT DATA SHEET

INFUSOMAT SP. LINE,SAFESET,PVC,LL,300CM (8270358SP)

Attribute Name	Attribute Value
Brand name	Infusomat®
Product name (HC)	Infusomat® Space Line SafeSet
Product type	Administration set
Application	IV-Standard
Type of administration set	Pressure
Container connector	IV standard spike
Length	300.00 cm
Length from pump outlet	200 cm
Inner-ø tube	3.00 mm
Outer-ø tube	4.10 mm
Type of filter	Fluid filter AirStop
Position of filter	Drip chamber
Pore size (Fluid filter)	15 µm
drops/ml	20
Air ventilation	vented

Attribute Name	Attribute Value
Flow control device	Roller clamp
Injection access	-
with AirStop	Yes
with PrimeStop	Yes
Number of back check valves	-
Number of SafSite®	-
Number of needle-free valves	-
Connection	Luer-Lock/male
Light-protected	No
Color of light protection	-
Pump compatibility (Type)	Infusomat® Space, Infusomat® fm/fmS
Length of pump segment	100.00 mm
max. pressure	2.0 bar
Material	PVC
Catalogue attribute	IV-Standard
Type of packaging	Blister
Product (unpack.) contains PVC	Yes
Product (unpack.) cont. Latex	No
Product (unpack.) cont. DEHP	No
Packaging contains PVC	No
Packaging contains Latex	No
Packaging contains DEHP	No
Product contains PVC	Yes
Product contains Latex	No
Product contains DEHP	No
Product contains Nickel	No
Product contains BPA	No
Disposable / Reusable	disposable
Sterility	sterile
Type of Sterilization	EO
Type of Sterilizat. (allowed)	EO
Priming volume	21.300 ml

Attribute Name	Attribute Value
Length (Inch)	118"