


PRODUCT DATA SPECIFICATION SHEET

PRODUCT REF CODE	SE-MIC002 (SY-LHC01E)			
PRODUCT DESCRIPTION	Light Handle Cover for Deroyal handle			
COMPOSITION	LDPE			
PRODUCT SIZE	Opening Size	Sheat Size		
	32mm (+0.2mm)	90mm		
UNIT OF MEASURE	Double sterile soft blister pack		IS PRODUCT STERILE?	Yes
RE-USE	NO		PHTHALATE STATUS	PHTHALATE FREE
CE CLASSIFICATION	Is		LATEX STATUS	LATEX FREE
PACKAGING	Inner box: 38*28*25cm, 80pc/box, 320pc/ctn, 8.4kg/ctn Outer case: 58*40*52cm			
LABELLING	IN ACCORDANCE WITH ALL RELEVANT SPECIFICATIONS AND STANDARDS INCLUDING THE MEDICAL DEVICES DIRECTIVE, 93/42/EEC, ISO15223-1:2012 AND ISO13485:2016			
SHELF LIFE	5 YEARS			
STORAGE CONDITIONS	IN CLEAN DRY WAREHOUSE IN AMBIENT CONDITIONS AND PROTECTED FROM SUNLIGHT.			
HANDLING	HANDLE WITH CARE			
METHODS OF STERILISATION	SUITABLE FOR ETHYLENE OXIDE STERILISATION			

ISSUE NUMBER		DATE	18-Feb-2021
SIGNED	Kevin Jiang	ROLE	Product manager

Intended Purpose and Using instruction:

A Lamp Handle Cover is a plastic cover intended to avoid direct touch to the surgical lamp handle when doctor adjust the lamp position causing across contamination. It like a flexible sleeve set in the Lamp Handle. By holding the sleeve, the doctor can easily slide the handle cover over the grip handle and adjust the handle cover in a secure position. The doctor can easily disengage the handle cover after the surgical procedure by pulling it off..