


## **Product Data Sheet**

**Item Description:** Sanibel™ insert foam tips





Frequency response measured of an IP30 in a 711 coupler using the Sanibel  $^{\rm TM}$  insert foam tip

Adult , medium foam tip	Pediatric, small foam tip
Ø approximately 13,7 mm	Ø approximately 9,7 mm
Packaging: Bag/50pcs - 8527918 Bag/100pcs - 8527920	Packaging: Bag/50pcs - 8527912 Bag/100pcs - 8527913

**Use:** AUD/IMP/ABR. Compatible with all types of insert phones



Conforms to Medical Device Regulation



LatexFree



Single patient use only!