

PRODUCT DATA SHEET

DAFILON BLUE 6/0 (0.7) 45CM DS16 (C0932116)



Attribute Name	Attribute Value
Aesculap Product Groups	Sutures
Product name	DAFILON
Color of thread	blue
Size of Thread (USP - Others)	6/0
Size of Thread (EP)	0,7
Number of threads per piece	1
Number of needles per piece	1
Length of thread	45 cm
Threads per peel	45 cm
Number of needles per thread	1
Needle 1	DS16
Needle 2	-
Needle Code	DS16
Needle Trade Name	Standard needle
Type of packaging	Direct Dispense Packaging
Loop	-
Take Off	-
Aro	-

Attribute Name	Attribute Value
Pledget	-
Type of Pledget Material	-
Axial	-
Target market	B.Braun standard
Observations	-
Degradation profile	non-absorbable
Coating of Thread	Uncoated
Structure of Thread	Monofilament
Origin of thread	Synthetic
Thread Material	polyamide
Angle of Needle 1	3/8 circle
Body of Needle 1	reverse cutting needle
Point of Needle 1	-
Length of Needle 1	16.0 mm
Thickness of Needle 1	normal
Easyblack	-
Surgical Step	Closure
Sales unit	box of 36 pieces
Product contains BBP	No
Product contains DBP	No
Product contains DEHP	No
Product contains DIDP	No
Product contains DINP	No
Product contains DOP	No
Product contains Latex	No
Product is powdered	No
Product contains PVC	No
Product contains chrome	Yes
Product contains Nickel	Yes
Sterility	sterile
Type of Sterilizat. (allowed)	for single use, not re-sterilizable
Material	stainless steel

Attribute Name	Attribute Value
Creation date	2001-02-05
Material Number	C0932116
Disposable / Reusable	disposable