



Image Coming Soon

17315

C-Arm tape drape set VE:25

QUANTITY: **25**

CONTAINS LATEX: **No**

STERILE: **Yes**

STERILIZATION METHOD: **EO**

MATERIAL: **Polyethylene**

GTIN PRIMARY: **50748426123341**

GTIN SECONDARY: **00748426123346**

GTIN SALES: **50748426123341**

COMMODITY CODE: **3926909790**

SHELF LIFE: **5 Years**

DEHP/PHTALATES: **No**

GMDN CODE: **12535**

CE NOTIFYING BODY: **CE0044**

MEDICAL CLASSIFICATION: **Is**

BOX LENGTH (MM): **460**

BOX WIDTH (MM): **295**

BOX HEIGHT: **275**

BOX M3: **0.0373175**

BOX WEIGHT (GRAM): **6000**

LABEL LANGUAGES: **BG,CS,DA,DE,EL,EN,ES,ET,FI,FR,HR,HU,IS,IT,LT,LV,NL,NO,PL,PT,RO,RU,SK,SL,SR,SV,TR**

2X STERILISATION ALLOWED?: **No**

ADDRESS OF LEGAL MANUFACTURER: **1 Ecolab Place, St. Paul, MN 55102 U.S.A.**

IFU: **No**

SINGLE USE: **Yes**

CASES PER PALLET CONVERSION: **20**

INFORMATION LAST UPDATED ON: **21 Jun 2022**
