



PC1295EU

Surgical / Intraoperative probe
cover

SIZE: **10 x 244 cm**

QUANTITY: **20**

CONTAINS LATEX: **No**

STERILE: **Yes**

STERILIZATION METHOD: **EO**

MATERIAL: **Ethyl Methyl Acrylate/anti-stat blend**

FOLDED/ROLLED: **Folded**

GTIN PRIMARY: **50748426099509**

GTIN SECONDARY: **00748426099504**

GTIN SALES: **50748426099509**

COMMODITY CODE: **3926909790**

INTENDED USE: **PC**

SHELF LIFE: **5 Years**

DEHP/PHTALATES: **No**

GMDN CODE: **44714**

CE NOTIFYING BODY: **CE0482**

MEDICAL CLASSIFICATION: **IIa**

BOX LENGTH (MM): **280**

BOX WIDTH (MM): **225**

BOX HEIGHT: **140**

BOX M3: **0.00882**

BOX WEIGHT (GRAM): **1600**

LABEL LANGUAGES: **BG,CZ,DE,DK,EE,EN,ES,FI,FR,GR,HR,HU,IS,IT,LT,LV,NL,NO,PL,PT,RO,RU,SE,SI,SK,SR,TR**

NON STERILE ALTERNATIVE: **Check with EU Marketing**

2X STERILISATION ALLOWED?: **No**

ADDRESS OF LEGAL MANUFACTURER: **Sorbonne Centre F20, Mosta, Technopark, Mosta, MST3000 Malta**

BISPHENOL A: **No**

IFU: **Yes**

SINGLE USE: **Yes**

CASES PER PALLET CONVERSION: **66**

INFORMATION LAST UPDATED ON: **21 Jun 2022**
