



## PC3787EU

Surgical / Intraoperative probe  
cover Straight tip

SIZE: **12 x 244 cm**

---

QUANTITY: **20**

---

CONTAINS LATEX: **No**

---

STERILE: **Yes**

---

STERILIZATION METHOD: **EO**

---

MATERIAL: **Ethyl Methyl Acrylate/anti-stat blend & PIP**

---

FOLDED/ROLLED: **Folded**

---

GTIN PRIMARY: **50748426099608**

---

GTIN SECONDARY: **00748426099603**

---

GTIN SALES: **50748426099608**

---

COMMODITY CODE: **3926909790**

---

INTENDED USE: **PC**

---

SHELF LIFE: **Expiration date must match the exp date on the 63132 Gel**

---

DEHP/PHTALATES: **No**

---

GMDN CODE: **44714**

---

CE NOTIFYING BODY: **CE0482**

---

MEDICAL CLASSIFICATION: **Ila**

---

BOX LENGTH (MM): **286**

---

BOX WIDTH (MM): **223**

---

BOX HEIGHT: **140**

---

BOX M3: **0.00892892**

---

BOX WEIGHT (GRAM): **1680**

---

COUNTRY OF ORIGIN: **Malta**

---

LABEL LANGUAGES: **BG,CZ,DE,DK,EE,EN,ES,FI,FR,GR,HR,HU,IS,IT,LT,LV,NL,NO,PL,PT,RO,RU,SE,SI,SK,SR,TR**

---

NON STERILE ALTERNATIVE: **PC3787N**

---

2X STERILISATION ALLOWED?: **No**

---

ADDRESS OF LEGAL MANUFACTURER: **Sorbonne Centre F20, Mosta, Technopark, Mosta, MST3000 Malta**

---

IFU: **Yes**

---

SINGLE USE: **Yes**

---

CASES PER PALLET CONVERSION: **156**

---

INFORMATION LAST UPDATED ON: **21 Jun 2022**

---