

**86715****Ultracover 44/300 mm VE: 100**SIZE: **44 x 300 mm**QUANTITY: **100**CONTAINS LATEX: **Yes**STERILE: **No**STERILIZATION METHOD: **Non Sterile**MATERIAL: **Latex**FOLDED/ROLLED: **Rolled**GTIN PRIMARY: **50748426094429**GTIN SALES: **50748426094429**COMMODITY CODE: **4014100000**INTENDED USE: **to cover an ultrasound transducer, used in natural body orifices, less than 1 hour**SHELF LIFE: **3 Years**DEHP/PHTALATES: **No**GMDN CODE: **44713**CE NOTIFYING BODY: **CE0482**MEDICAL CLASSIFICATION: **Ila**BOX LENGTH (MM): **290**BOX WIDTH (MM): **200**BOX HEIGHT: **125**

BOX M3: **0.00725**

BOX WEIGHT (GRAM): **1150**

LABEL LANGUAGES: **DA,DE,ES,FI,FR,IT,NL,PL,PT,SE**

NON STERILE ALTERNATIVE: **0**

2X STERILISATION ALLOWED?: **No**

ADDRESS OF LEGAL MANUFACTURER: **Hekkehorst 24, 7207 BN Zutphen, The Netherlands**

IFU: **Yes**

SINGLE USE: **Yes**

SPECIAL STORAGE CONDITIONS?: **Keep Away From Sunlight Symbol**

CASES PER PALLET CONVERSION: **240**

INFORMATION LAST UPDATED ON: **21 Jun 2022**
