



# 17558

## Light handle cover, large

SIZE: **Large**

---

QUANTITY: **200**

---

CONTAINS LATEX: **No**

---

STERILE: **Yes**

---

STERILIZATION METHOD: **EO**

---

MATERIAL: **Polyethylene tubing, Polyethylene ring**

---

GTIN PRIMARY: **50748426084420**

---

GTIN SECONDARY: **00748426084425**

---

GTIN SALES: **50748426084420**

---

COMMODITY CODE: **3926909790**

---

INTENDED USE: **covering equipment**

---

SHELF LIFE: **5 Years**

---

DEHP/PHTHALATES: **No**

---

GMDN CODE: **44977**

---

CE NOTIFYING BODY: **CE0482**

---

MEDICAL CLASSIFICATION: **Is**

---

BOX LENGTH (MM): **600**

---

BOX WIDTH (MM): **400**

---

BOX HEIGHT: **300**

---

BOX M3: **0.072**

---

BOX WEIGHT (GRAM): **7850**

---

LABEL LANGUAGES: **BG,CZ,DE,DK,ES,FI,FR,GR,HR,HU,IT,NL,NO,PL,PT,RU,SE,TR**

---

NON STERILE ALTERNATIVE: **3317558**

---

2X STERILISATION ALLOWED?: **No**

---

ADDRESS OF LEGAL MANUFACTURER: **Hekkehorst 24, 7207 BN Zutphen, The Netherlands**

---

BISPHENOL A: **No**

---

IFU: **No**

---

PVC: **No**

---

REACH COMPLIANCE: **Yes**

---

SINGLE USE: **Yes**

---

CASES PER PALLET CONVERSION: **18**

---

INFORMATION LAST UPDATED ON: **21 Jun 2022**

---