

Article no. 70020
Category Cellulose wadding
Product Cellulose wadding 40 x 60 cm, bleached
End-use Absorbing liquids, wiping off surfaces, etc.



1. Ideal and hygienic to use for absorbing liquids, wiping off surfaces, etc.
2. With a wide usage range e.g. health care, dental, veterinary, etc.
3. With high absorption capacity.
4. Supplied in plain sheets.
5. Made of high-quality Scandinavian tissue, bleached.
6. 100% virgin fibers of TCF-pulp (Totally Chlorine Free).
7. Certified according to Nordic Swan ecolabel.
8. 5 kg/carton.
9. CE-marked product. Fulfils the requirements of (EU) MDR 2017/745.
10. Packed in cartons on pallets.
11. The pallets are stretched-wrapped with transparent PE-film (polyethylene).
12. The cartons are made of cardboard and could be incinerated or recycled under standard conditions as cardboard.
13. The stretch film could be incinerated during standard conditions, or recycled as standard plastic.
14. Produced in Sweden.



Size (cm)	Material	Units/ inner package Inner packages/ transport unit	Transport units/ pallet Height of pallet (m)	Gross/Net weight carton (kg)	Measure of carton LxWxH (cm)	Pallet (EUR) size (mm)
40 x 60	Tissue	1 1	24 1,020	5,0 4,6	60 42 14,5	1200 x 800

*Finess Hygiene is certified in accordance with ISO 9001:2015, ISO 14001:2015 and ISO 13485:2016.
 Our CE-marked products fulfils the requirements of (EU) MDR 2017/745.*

2022-12-29 ver. 04