


Verifine®

Insulin Syringes

**YOUR IDEAL CHOICE FOR
INSULIN SYRIGES !**



FEATURES AND BENEFITS



Comfortable

- Ultra thin, sharp and specially lubricated needles

Readable

- Graduation lines in bold print for easy readability

Reliable

- Quality design to avoid waste of insulin
- Quality design for dosing precision

AVAILABLE IN INDIVIDUAL BAGS AND BAGS OF MORE THAN ONE PIECE



Type 7

Individual bags



Type 8

Bags of more than one piece
The number of pieces can be customized

ORDERING INFORMATION

Insulin Syringe	REF Type 7	REF Type 8	Needle Size	Syringe size	Qty/Box	Qty/Case
U-40	ISY-7053008-40	ISY-8053008-40	30G x 5/16"	0.5ml	100	2000
	ISY-7053012-40	ISY-8053012-40	30G x 1/2"	0.5ml	100	2000
	ISY-7052912-40	ISY-8052912-40	29G x 1/2"	0.5ml	100	2000
	ISY-7103108-40	ISY-8103108-40	31G x 5/16"	1ml	100	2000
	ISY-7103112-40	ISY-8103112-40	31G x 1/2"	1ml	100	2000
	ISY-7103008-40	ISY-8103008-40	30G x 5/16"	1ml	100	2000
	ISY-7103012-40	ISY-8103012-40	30G x 1/2"	1ml	100	2000
	ISY-7102912-40	ISY-8102912-40	29G x 1/2"	1ml	100	2000
U-100	ISY-7033108-100	ISY-8033108-100	31G x 5/16"	0.3ml	100	2000
	ISY-7033008-100	ISY-8033008-100	30G x 5/16"	0.3ml	100	2000
	ISY-7032912-100	ISY-8032912-100	29G x 1/2"	0.3ml	100	2000
	ISY-7053105-100	ISY-8053105-100	31G x 3/16"	0.5ml	100	2000
	ISY-7053108-100	ISY-8053108-100	31G x 5/16"	0.5ml	100	2000
	ISY-7053008-100	ISY-8053008-100	30G x 5/16"	0.5ml	100	2000
	ISY-7053012-100	ISY-8053012-100	30G x 1/2"	0.5ml	100	2000
	ISY-7052912-100	ISY-8052912-100	29G x 1/2"	0.5ml	100	2000
	ISY-7103105-100	ISY-8103105-100	31G x 3/16"	1ml	100	2000
	ISY-7103108-100	ISY-8103108-100	31G x 5/16"	1ml	100	2000
	ISY-7103008-100	ISY-8103008-100	30G x 5/16"	1ml	100	2000
	ISY-7103112-100	ISY-8103112-100	31G x 1/2"	1ml	100	2000
	ISY-7103012-100	ISY-8103012-100	30G x 1/2"	1ml	100	2000
	ISY-7102912-100	ISY-8102912-100	29G x 1/2"	1ml	100	2000

✘ If you have further questions or need other specifications, please contact us.