

# DANATUBE

Elastisk tubegaze,  
lavet af bambus.



LANGVARENDE FRISKHED




VIDUNDERLIG BLØD



FREMRAGENDE VARMEREGULERING



MILJØVENLIG



MEDIPLAST kan tilbyde DANATUBE, en elastisk tubeforbinding med god luftgennemtrængelighed.

## DANATUBE BAMBUS TUBEGAZE

LÆNGDE: 20 METER USTRAKT

Æske Trsp.karton

6050731	Size 0,1 = 1,5 cm, fingers and toes	1	64
6050732	Size 12 = 2,5 cm, fingers and toes with cotton wool wadding	1	33
6050733	Size 34 = 4 cm, hand, small arms and childrens legs	1	29
6050734	Size 56 = 6 cm, large arms and normal legs	1	60
6050735	Size 78 = 8 cm, head, knee and small thighs	1	48
6050736	Size T1 = 10 cm, trunk bandage for children and large head	1	39
6050738	Size = 12 cm, trunk bandage -medium	1	36
6050737	Size T2 = 16 cm, trunk bandage for adult	1	24

Æske = mindste kvantitet/bestilling

OEKO-TEX®  
CONFIDENCE IN TEXTILES  
STANDARD 100



FREMSTILLET I DANMARK

Indikation: Til dermatologisk brug. Under gips og skinner eller til fiksering af sårforbindingsbandage. Må kun anvendes på ubeskadiget hud.