



Description

Reference

- 0121-12 MIC-KEY* Extension set with right-angle connector – 12 inch - 30 cm
- 0121-24 MIC-KEY* Extension set with right-angle connector – 24 inch – 60 cm
- 0122-02 MIC-KEY* Medication set with right-angle connector – 2 inch – 5 cm
- 0123-12 MIC-KEY* Bolus Extension set with straight connector – 12 inch – 30 cm
- 0123-24 MIC-KEY* Bolus Extension set with straight connector – 24 inch – 60 cm
- 0124-12 MIC-KEY* Extension set Single Ended Funnel, with right-angle connector – 12 inch – 30 cm
- 0124-24 MIC-KEY* Extension set Single Ended Funnel, with right-angle connector – 24 inch – 60 cm
- 0126-12 MIC-KEY* Extension set Single Ended Luer Lock with right-angle connector – 12 inch – 30 cm
- 0126-24 MIC-KEY* Extension set Single Ended Luer Lock with right-angle connector – 24 inch – 60 cm

The MIC-KEY* Extension Sets are intended to connect to the MIC-KEY* Low-Profile Feeding Tubes, Luer Slip or Catheter tip syringes and feed administration bags.

The MIC-KEY* Extension set has on the patient end the SECUR-LOK* connector and on the other side either a funnel or Luer slip connector.

The MIC-KEY* Extension sets (except the Medication set) have a Quick Release Clamp.

Indication

For the administration of enteral nutrition and/or medication only.

Single patient use only.

The MIC-KEY* Extension sets are disposables. For hygienic and functionality reasons, it is recommended to replace the MIC-KEY* Extension Sets at least every 2 weeks.

Counter indication

Consult the direction for use for additional information and warnings.

Main materials

Extension set tubing and connector: PVC with adipate as softener.

Patient connector: ABS. Polyurethane

Clamp: Polypropylene

Not formulated with Natural Rubber Latex.

Not formulated with DEHP / DBP / BBP / Bisphenol A.



REF 0121-12

REF 0121-24



REF 0122-02



REF 0123-12

REF 0123-24



REF 0124-12

REF 0124-24



REF 0126-12

REF 0126-24

Packaging

Each Extension set is packed in a plastic pouch.

5 plastic pouches are placed in a carton box (shipping case).

Box size: (LxWxD) 15,9 x 15,6 x 5,0 cm

Box weight: 172g (0,38lbs)

Bar coding: GTIN symbology, linear, on unit packaging and carton box.

Item #	Single Unit	Carton box
0121-12	00350770949277	10350770949274
0121-24	00350770949291	10350770949298
0122-02	00350770949321	10350770949328
0123-12	00350770949338	10350770949335
0123-24	00350770949352	10350770949359
0124-12	00350770949376	10350770949373
0124-24	00350770949390	10350770949397
0126-12	00350770949437	10350770949434
0126-24	00350770949444	10350770949441

*Registered trademark or trademark of Avanos Medical Inc.

The dimensions and properties listed above can vary within pre-established specifications. This document was created using the most recent information. In the interest of continuous improvement, the characteristics of the product may change without prior notice.



Manufacturing

Manufactured in Mexico.
The manufacturing site's quality is ISO13485 compliant.

Regulatory information

Product CE marked as per 93/42/EEC Directive.
Class of the device: IIa.
Notified body: BSI (2797).

Sterilisation

These products are non-sterile.

Storage

Store in a dry and cool place, away from sources of heat and radiation.
Keep as much as practicably possible in its carton box.

Shelf life

5 years, from date of manufacturing.