

**3823****Banded bag, circular sewn elastic**SIZE: **191 x 94 cm**QUANTITY: **20**CONTAINS LATEX: **No**STERILE: **Yes**STERILIZATION METHOD: **GAMMA**MATERIAL: **Polyethylene**GTIN PRIMARY: **50748426010757**GTIN SECONDARY: **00748426010752**GTIN SALES: **50748426010757**COMMODITY CODE: **39269097**SHELF LIFE: **4 Years**DEHP/PHTALATES: **No**GMDN CODE: **12535**CE NOTIFYING BODY: **CE0044**MEDICAL CLASSIFICATION: **Is**BOX LENGTH (MM): **540**BOX WIDTH (MM): **440**BOX HEIGHT: **290**BOX M3: **0.068904**

BOX WEIGHT (GRAM): **5750**

---

LABEL LANGUAGES: **CS,DA,DE,EL,EN,ES,FI,FR,HU,IT,NL,PL,PT,SV**

---

2X STERILISATION ALLOWED?: **No**

---

ADDRESS OF LEGAL MANUFACTURER: **1 Ecolab Place, St. Paul, MN 55102 U.S.A.**

---

BISPHENOL A: **No**

---

IFU: **Yes**

---

REACH COMPLIANCE: **Yes**

---

SINGLE USE: **Yes**

---

CASES PER PALLET CONVERSION: **48**

---

INFORMATION LAST UPDATED ON: **21 Jun 2022**

---