



PC1293

Surgical / Intraoperative probe cover

SIZE: **8 x 244 cm**

QUANTITY: **20**

CONTAINS LATEX: **No**

STERILE: **Yes**

STERILIZATION METHOD: **EO**

MATERIAL: **Ethyl Methyl Acrylate/anti-stat blend**

FOLDED/ROLLED: **Folded**

GTIN PRIMARY: **50748426047876**

GTIN SECONDARY: **00748426047871**

GTIN SALES: **50748426047876**

COMMODITY CODE: **3926909790**

INTENDED USE: **A device intended to be used as a physical barrier for protection against the effects of environmental exposure (e.g., body fluids, gels) and/or to maintain the required hygienic level of various dia**

SHELF LIFE: **Expiration date must match the exp date on the 63132 Gel**

DEHP/PHTALATES: **No**

GMDN CODE: **44714**

CE NOTIFYING BODY: **CE0044**

MEDICAL CLASSIFICATION: **Ila**

BOX LENGTH (MM): **290**

BOX WIDTH (MM): **230**

BOX HEIGHT: **140**

BOX M3: **0.009338**

BOX WEIGHT (GRAM): **1500**

COUNTRY OF ORIGIN: **Dominican Republic**

LABEL LANGUAGES: **CZ,DE,DK,EN,ES,FR,GR,HU,IT,NL,PL,PT,RU,SE**

NON STERILE ALTERNATIVE: **Check with EU Marketing**

2X STERILISATION ALLOWED?: **Yes**

ADDRESS OF LEGAL MANUFACTURER: **1 Ecolab Place, St. Paul, MN 55102 U.S.A.**

BISPHENOL A: **No**

IFU: **Yes**

PVC: **No**

REACH COMPLIANCE: **Yes**

SINGLE USE: **Yes**

CASES PER PALLET CONVERSION: **156**

INFORMATION LAST UPDATED ON: **21 Jun 2022**
