



PC1230EU

***GENERAL PURPOSE PROBE
COVER, 3 X 24**

SIZE: **3 x 24"**

QUANTITY: **100**

CONTAINS LATEX: **No**

STERILE: **Yes**

STERILIZATION METHOD: **EO**

MATERIAL: **Ethyl Methyl Acrylate/anti-stat blend**

FOLDED/ROLLED: **Folded**

GTIN PRIMARY: **50748426088084**

GTIN SECONDARY: **00748426088089**

GTIN SALES: **50748426088084**

COMMODITY CODE: **3926909790**

INTENDED USE: **PC**

SHELF LIFE: **5 Years**

DEHP/PHTALATES: **No**

GMDN CODE: **44713**

CE NOTIFYING BODY: **CE0482**

MEDICAL CLASSIFICATION: **Ila**

BOX LENGTH (MM): **460**

BOX WIDTH (MM): **300**

BOX HEIGHT: **150**

BOX M3: **0.0207**

BOX WEIGHT (GRAM): **1750**

LABEL LANGUAGES: **BG,CZ,DE,DK,EE,EN,ES,FI,FR,GR,HR,HU,IS,IT,LT,LV,NL,NO,PL,PT,RO,RU,SE,SI,SK,SR,TR**

NON STERILE ALTERNATIVE: **0**

2X STERILISATION ALLOWED?: **Yes**

ADDRESS OF LEGAL MANUFACTURER: **Sorbonne Centre F20, Mosta, Technopark, Mosta, MST3000 Malta**

IFU: **Yes**

SINGLE USE: **Yes**

CASES PER PALLET CONVERSION: **156**

INFORMATION LAST UPDATED ON: **21 Jun 2022**
