

## Technical Sheet

### PRO-Breathe Nasal Cannula

#### Key Features

- Over the ear design improves patient comfort during long term oxygen use
- Sliding toggle adjusts under the chin for best fit
- Soft, comfortable and flexible nasal prongs in a choice of styles
- Transparent with crush/kink resistant tubing for maximum patient safety
- Single use to help to prevent cross infection and contamination
- 

#### Intended Purpose

To be inserted into the nostrils of a patient to administer oxygen (O2).

Stock Code	Style	Nominal Size	Nominal Length (O2 Tubing)
PBNC10121	Straight Prongs	Adult	2.1m
PBNC10221	Straight Prongs	Paediatric	2.1m
PBNC20121	Curved Prongs	Adult	2.1m
PBNC20221	Curved Prongs	Paediatric	2.1m
PBNC20321	Curved Prongs	Infant	2.1m
PBNC20421	Curved Prongs	Neonate	2.1m
PBNC20150	Curved Prongs	Adult	0.5m

Material	Phthalates	Latex	Sterility	Shelf Life	Reusability	Use By
PVC	Yes	No	Non-Sterile	5 Years	Single Use	Yes

Device Class	Class IIa	Device Classification	Rule 2, MDD 93/42/EEC
Manufacturer	Well Lead Medical Co.Ltd	Authorised Representative	Shanghai International Holding Corp.
Importer Name	PROACT Medical Ltd	Country of Origin	China
Intrastat Code	9018390000	GMDN Reference	35201