




CODAN Spikes

CODAN Spike



1. A compact and user friendly system for preparing and withdrawing medicines. Specially shaped piercing device. The pressure-resistant hydrophobic 0,2 µm filter guarantees sterile ventilation and reliably excludes toxic aerosols during the preparation of medicines. The high-quality housing materials ensures that there is no loss of active substances during the preparation and withdrawal of medicines. The large surface area of the filter membrane ensures easy fluid transfer without high working pressure. Furthermore, fluids can be withdrawn without leaving significant residues.



Disclaimer:

The above is indicated only for information. This document is not a specification. The values given are subject to the usual fluctuations. Please follow the instructions for use by the medication manufacturer and the actual SPC (Summary of Product Characteristics). Common hygiene measures and working methods are to be followed at all times.

Subject to alterations without prior notice!

REF	Product description	Packaging	GTIN	Filling/Storage volume
16.5260	Withdrawal spike CODAN Spike, 0.2 µm air filter	1 50 600 14400	04055698001972 04055698001989 04055698001996 -	



CODAN Spikes

CODAN Spike

REF	Product description	Packaging	GTIN	Filling/Storage volume
16.5274	Withdrawal spike 5x CODAN Spike, 0.2 µm air filter	5	04055698009428	
		250	04055698009435	
		18000	-	
16.5276	Withdrawal spike 25x CODAN Spike, 0.2 µm air filter	25	04055698015160	
		500	04055698015177	
		36000	-	