

InsulinSaver – Smart Insulin Thermometer



InsulinSaver is designed to be a clear and user-friendly aid. InsulinSaver continuously analyzes the temperature that the insulin has been exposed to. It flashes red if the temperature is or has been above +30°C, and it flashes blue for a temperature below +2°C. InsulinSaver is light and easy to carry along with your insulin. The product is designed so that it easily fits it into, for example, your injection case or by your insulin packaging.

- InsulinSaver is very easy to use, “set & forget”
- Use InsulinSaver for increased safety when traveling
- Forget the worry when storing your insulin in the refrigerator

Specifications:

Accuracy	± 0,5 °C
Verification method of accuracy	Measuring equipment calibrated according to IPTS-68
Temperature sensor	Negative Temperature Coefficient Thermistor
Upper temperature limit	+30 °C
Lower temperature limit	+2 °C
Battery type	Lithium Manganese Dioxide
Alarm type	LED
Waterproof	YES
IP class	IP68
Possible to restore	YES
Battery time	12-18 months <i>18 months continuous operation according to test procedure</i>
Dimensions	5,5cm x 2,5cm x 0,4cm
EAN / PN	7300009080069