


The Bag - the safest way to store your insulin at the perfect temperature.

Compared to..

.. other insulin bags on the market, our bag is the **best in test**. It maintains the right temperature longer than any other bag – regardless of humidity. And yes, we've tested it more than once.



The best thing?

You don't have to freeze it, soak it in water or charge it. The Bag works all by itself. No stress and no extra planning - **just grab it and go.**

How does it work?

The Bag contains a special cooling gel material that activates at temperatures above 25°C/77°F, transitioning from a solid to liquid state as it absorbs heat. Keeping your insulin below 30°C/86°F - even in hot conditions.

Specifications

Recharging method

Cooling naturally below 24°C/75°F or refrigerating (optional)

Max temperature protection

54°C/129°F

Duration at 35°C/95°F

8 hours +

Self-cooling insulin bag.



Insulin safety made easy.

